



Decal Stat™ Decalcifier

Technical Data Sheet & IFU

Decal Stat™ is a hydrochloric acid solution which contains a chelating agent. Decal Stat™ is formulated to work faster than Original Decal® product. Can be used for the rapid decalcification of all types of bone and it works well with routine and all types of special stains. Cellular and nuclear details remain crisp and clear when the procedure is carefully monitored. Minimal distortion of tissue and cellular structure after the decalcification process make this an ideal product for your calcified tissues when Immunohistochemistry (IHC) procedures are not being considered. If IHC is part of the diagnostic workup, then Formalin 4™, Formalin 2000™, or Immunocal™ should be considered.

APPLICATIONS/INTENDED USE

Decal Stat™ may be used to process all types of calcified histological specimens. Formulated to work faster than original Decal®. All types of routine and special stains work flawlessly as well.

Fix tissue in 10% Neutral Buffered Formalin or other suitable tissue fixative such as Bouins, B5, or Zinc Formalin. The volume of fixative solution should be 20 times the tissue volume.

RESULTS

To determine the end point of the decalcification process use either x-ray or chemical end-point determination techniques. Do not use probes, needles, scalpels or bending as this may cause physical damage to the specimen.

PROCEDURE

1. Specimens must be fixed before exposure to an acid solution.
2. It is best to suspend the specimen so it is not in contact with any of the surfaces of the container. This allows exposure to all specimen surfaces and allows the precipitated calcium salts to sink to the bottom of the container.
3. Small specimens should not be left in the solution overnight. If decalcification process is incomplete, wash it in water and return to fixative. Wash again before returning to the decalcification solution to complete the process.
4. Continue Step #3 until specimen is completely decalcified.
5. If decalcification time exceeds 24 hours, it is best to replace the decalcifying solution with fresh solution daily.
6. **Please Note:** If you plan to stain the section with a Potassium Ferrocyanide / HCL stain (Iron stain), a minimum 10 minute rinse after decalcification is recommended.
7. Process tissues, embed and cut per Lab SOP (Standard Operating Procedure)

STORAGE/ SAFETY

Storage: Room Temperature
Refer to SDS for details-

Caution:

Specimens must be washed in running water before and after exposure to acid solutions, especially hydrochloric acid. The combination of formalin and hydrochloric acid can create the formation of **bis-chloromethyl ether which is a known carcinogen.**

ORDERING INFORMATION

Product		Product Number
Decal Stat™	32 ounces	D1212-32
Decal Stat™	1 gallon	D1212-1
Decal Stat™	5 gallon	D1212-5

TECHNICAL INFORMATION

1. If decalcification process was incomplete, surface decalcification techniques may be used.
2. To surface decal: Place a small dish of Decal Stat™ on ice. Place the face of the block in the dish for 5-10 minutes. Rinse the block in cold water. Icing tends to make the block harder and the water shed tends to soften the tissue face. Icing will greatly reduce the amount of chattering, especially in large blocks. Be certain to rinse the block before it is placed on the microtome (decal solution is corrosive)
3. Remember that bone has all of the same tissue elements as any other tissue and in order to demonstrate them effectively the specimens must be well fixed and carefully monitored during the decalcification process.

CONTACT INFORMATION

StatLab Medical Products, Inc.
2090 Commerce Drive
McKinney, Texas 75069
Phone: 800-442-3573 Fax: 972-436-1369
www.statlab.com

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or give credit for the purchase price. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for their own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any violation of any law or government regulation.

21-May-2015