



American MasterTech
scientific laboratory supplies

Carnoy's Fluid

Ready to Use Histological Tissue Fixative

(Rev. 3/8/18)

Carnoy's Fluid The product and its purpose:

Carnoy's Fluid is a Histological tissue fixative routinely used preserving RNA, nissl granules, and glycogen. Carnoy's Fluid is intended for In Vitro Diagnostic use. Tissue fixed with Carnoy's Fluid can be used for routine Special Stains and Immunohistochemical (IHC) stains. Carnoy's Fluid is compatible with routine processing alcohols. Active ingredients include; Ethanol (<65%), Chloroform (<30%), and Acetic Acid (<15%).

Carnoy's Fluid, Instructions for use:

Use Carnoy's Fluid as you would any alcohol based fixative. Most surgical tissue specimens require 6 to 24 hours for complete fixation. Large tissues may take longer to ensure complete penetration and fixation.

Carnoy's Fluid, Storage and Stability:

Carnoy's Fluid should be stored at 15-30 C and is stable under normal laboratory conditions. See product label for shelf life.

Carnoy's Fluid, Precautions:

See Safety Data Sheet (SDS) for warnings and precautions.

Carnoy's Fluid, Disposal:

Dispose Carnoy's Fluid in accordance with all local, state, federal, environmental and/or EU regulations.

CE MARKINGS AND DESIGNATIONS:

| | | | | | | |
|------------|------------------------------------|-----------|-----------------------------------|---------------|----------------|--|
| REF | Catalogue Number | | Temperature Limitation | | Manufacturer | American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136 |
| LOT | Batch Code | | Use By | EC REP | Representative | Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands |
| IVD | In Vitro Diagnostic Medical Device | | Consult Instructions Prior to Use | GHS02 | Flammable | CE |
| GHS06 | Toxic | GHS07 | Irritant | GHS08 | Health Hazard | |

Copyright © 2015-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.