



American MasterTech  
scientific laboratory supplies

# Frozen Section Fixative

## Ready to Use Histological Tissue Fixative

(Rev. 3/8/18)

### Frozen Section Fixative, The product and its purpose:

Frozen Section Fixative is a Histological tissue fixative intended for In Vitro Diagnostic use. Frozen tissue sections fixed with Frozen Section Fixative can be used for routine Special Stains and Immunohistochemical (IHC) stains. Antigen retrieval will be needed for IHC staining. Active ingredients include; Isopropyl Alcohol (<80%), Glyoxal (<10), and Ethanol (<5%).

### Frozen Section Fixative, Instructions for use:

After cutting the frozen section and mounting on a slide, place the slide in a coplin jar containing Frozen Section Fixative for 15-30 seconds. Next, carefully dip the slide 10 times in a coplin jar of water and stain as usual.

### Frozen Section Fixative, Storage and Stability:

Frozen Section Fixative should be stored at 15-30 C and is stable under normal laboratory conditions. See product label for shelf life.













### Frozen Section Fixative, Precautions:

See Safety Data Sheet (SDS) for warnings and precautions.

### Frozen Section Fixative, Disposal:

Dispose of or neutralize Frozen Section Fixative in accordance with all local, state, federal, environmental and/or EU regulations.

### CE MARKINGS AND DESIGNATIONS:

	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
	Batch Code		Use By		Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS02	Flammable	
 GHS07	Irritant	 GHS08	Health Hazard			

Copyright © 2015-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.