



## CD3 (LN10)

### Mouse Monoclonal Antibody

MM150-05

MM150-10

MM150-25

MM150-100

Document #: IFU-MM150-CD3\_LN10

Release Date: 12/01/17, IFU-006 Rev A

Immunogen	Clone	Species	Isotype	Primary Antibody Diluent
Prokaryotic recombinant protein corresponding to C-terminal region of CD3.	LN10	Mouse	IgG1	ACR-018

Lot Specific Ig concentration available upon request.

Catalog #	Description
MM150-05	0.5 mL concentrated antibody for use with StatLab Medical Ultra High Def Polymer Detection System on manual or automated IHC Systems
MM150-10	10ml barcoded ready to use antibody for use with StatLab Medical Ultra High Def Polymer Detection System on the automated IHC System
MM150-25	25ml ready to use antibody for use with StatLab Medical Ultra High Def Polymer Detection System on the automated IHC System
MM150-100	100ml ready to use antibody for use with StatLab Medical Ultra High Def Polymer Detection System on the automated IHC System

#### Intended Use

For In Vitro Diagnostic Use. This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. Clinical interpretation of staining results should be accompanied by histological studies with proper controls. Patients' clinical histories and other relevant diagnostic tests should be utilized by a qualified person(s) when evaluating and interpreting results.

#### Summary and Explanation

The Clone LN10 is specific for the non-glycosylated epsilon chain of the human CD3 molecule. Clone LN10 recognizes T cells in thymus, bone marrow, peripheral lymphoid tissue and blood and is a pan T cell marker. The CD3 molecule consists of five different polypeptide chains with molecular weights ranging from 16 to 28 kD. The CD3 antigen is first detected in early thymocytes and its appearance probably represents one of the earliest signs of commitment to the T cell lineage.

#### Format

This product is supplied as a tissue culture supernatant and contains sodium azide as a preservative.

#### Principles of the Procedures

Antigen detection by immunohistochemistry (IHC) is a two-step process involving first, the binding of a primary antibody to the antigen of interest, and second, the detection of bound antibody by a chromogen. The primary antibody may be used in IHC using manual techniques or using Automated Staining System.

#### Dilution of Primary Antibody

StatLab Medical Products ready to use antibodies have been optimized for use with the recommended StatLab Medical Products Detection System and do not require further dilution. Further dilution may result in loss of sensitivity. The user must validate any such change.

#### Materials Required But Not Provided

Some of the reagents and materials required for IHC are not provided. Pretreatment reagents, detection systems, control reagents and other ancillary reagents are available from StatLab Medical Products. Please refer to the StatLab Medical Products website at: [www.statlab.com](http://www.statlab.com)

#### Storage and Handling

Store at 2-8°C. This antibody is suitable for use until the expiration date when stored at 2-8°C. Do not use product after the expiration date printed on vial. If reagents are stored under conditions other than those specified here, they must be verified by the user. Diluted reagents should be used promptly. Unused portions of antibody preparation should be discarded after one day.

The presence of precipitate or an unusual odor indicates that the antibody is deteriorating and should not be used.

Positive and negative controls should be run simultaneously with all patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact StatLab IHC Technical Support via email at: [ihctech@statlab.com](mailto:ihctech@statlab.com) or call us at (800) 442-3573.

#### Specimen Collection and Preparation

Tissues fixed in 10% formalin are suitable for use prior to paraffin embedding. Consult references (Kiernan, 1981; Sheehan & Hrapchak, 1980) for further details on specimen preparation.

The user is advised to validate the use of the products with their tissue specimens prepared and handled in accordance with their laboratory practices.

#### Precautions

This antibody contains less than 0.1% sodium azide. Concentrations less than 0.1% are not reportable hazardous materials according to U.S. 29 CFR 1910.1200, OSHA Hazard Communication and EC Directive 91/155/EC. Sodium azide (NaN<sub>3</sub>) used as a preservative is toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with large volumes of water to prevent azide build-up in plumbing. (Center for Disease Control, 1976, National Institute of Occupational Safety and Health, 1976). Specimens, before and after fixation and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Microbial contamination of reagents may result in an increase in nonspecific staining. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change. The SDS is available upon request.

#### Treatment of Tissues Prior to Staining

Pretreatment of tissues, if any, should be done as suggested: Place the slides in the recommended Antigen Retrieval Solution using an appropriate retrieval/pressure cooker system. Set the temperature for 15-minute incubation at "High Pressure".

Formatted: Normal, Left, Space After: 0 pt



Manufactured by StatLab Medical Products  
2090 Commerce Drive  
McKinney, Texas 75069  
(800) 442-3573  
[www.StatLab.com](http://www.StatLab.com)

IVD



Allow slides to cool down for 20 minutes prior to staining.

Staining Procedure

Refer to the following table for conditions specifically recommended for this antibody. Refer to the StatLab Ultra High Def–Two Step Detection System for guidance on specific staining protocols or other requirements.

Parameter	StatLab Recommendations
Positive Control	Tonsil
Pretreatment	EDTA Buffer pH 8.0
Concentrated Dilution	1:50-1:80
Incubation Time & Temperature	30 min @ RT
Detection System	StatLab Medical Ultra High Def Polymer Detection System – Two Step Detection system
Tissue Type	FFPE

Quality Control

Refer to CLSI Quality Standards for Design and Implementation of Immunohistochemistry Assays; Approved Guideline-Second edition (I/LA28-A2) CLSI Wayne, PA USA (www.clsi.org). 2011.

Troubleshooting

StatLab Headquarters: 2090 Commerce Drive, McKinney, TX 75069.  
Email our team at [jhctech@statlab.com](mailto:jhctech@statlab.com)  
Call at (800) 442-3573.

Cellular Localization

Cell membrane

Limitations of the Procedure

IHC is a complex technique involving both histological and immunological detection methods. Tissue processing and handling prior to immunostaining can also cause inconsistent results. Variations in fixation and embedding or the inherent nature of the tissue may cause variations in results (Nadji and Morales, 1983). Endogenous peroxidase activity or pseudoperoxidase activity in erythrocytes and endogenous biotin may cause non-specific staining depending on detection system used. Tissues containing Hepatitis B surface Antigen (HBsAg) may give false positive with horseradish peroxidase systems (Omata et al, 1980). Improper counterstaining and mounting may compromise the interpretation of results.

Performance Characteristics

The optimum antibody dilution and protocols for a specific application can vary. These include, but are not limited to: fixation, heat-retrieval method, incubation times, and tissue section thickness and detection kit used. Due to the superior sensitivity of these unique reagents, the recommended incubation times and titers listed are not applicable to other detection systems, as results may vary. The data sheet recommendations and protocols are based on exclusive use of products manufactured for StatLab. Ultimately, it is the responsibility of the investigator to determine optimal conditions. These products are tools that can be used for interpretation of morphological findings in conjunction with other diagnostic tests and pertinent clinical data by a qualified pathologist.

References

I. Steward M, Bishop R, Piggott NH et al. Production and characterization of a new monoclonal antibody effective in recognizing the CD3 T-cell associated antigen in formalin-fixed embedded tissue. Histopathology 1997; 30, 16–22.

II. Williamson SLH, Steward M, Milton I et al. New monoclonal antibodies to the T cell antigens CD4 and CD8: Production and characterization in formalin-fixed, paraffin-embedded tissue. American Journal of Pathology 1998; 152 (6), 1421–1426.

III. Leong FJW-M and Leong AS-Y. Essential markers in malignant lymphoma: a diagnostic approach. The Journal of Histotechnology 2002; 25(4), 215–227.

IV. Krynitz B, Rozell B and Lindelof B. Differences in peritumoral inflammatory skin infiltrate between squamous cell carcinomas in organ transplant recipients and immunocompetent patients. Acta Dermato Venerologica. 2010; 90:379-385.