

Fontana-Masson Stain Kit (For Argentaffin Cells and Melanin)


Description: The Fontana-Masson Stain Kit is intended for use in the histological visualization of Argentaffin cells and Melanin in paraffin or frozen sections. In addition, the Fontana-Masson stain has been reported to be useful in identifying Capsule-Deficient Cryptococcus Neoformans and typical Cryptococcus Neoformans.

Uses/Limitations: Not to be taken internally. For In-Vitro Diagnostic use only. Histological applications. Do not use past expiration date. Use caution when handling reagents. Non-Sterile.

Control Tissue: Any paraffin embedded tissue that contains Hair Follicles or Skin for Melanin. Small Intestine for Argentaffin.

Results:

Argentaffin Cells:	Black
Melanin:	Black
Cryptococci Cell Wall:	Black
Nuclei:	Red
Cytoplasm:	Light Pink



Kit Contents:

<u>Item #</u>	<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
SSC-GCB125	Gold Chloride Solution (0.2%)	125 ml	2-8°C
SSC-SNX009	Silver Nitrate Solution (10%)	5 x 9ml	2-8°C
SSC-STB125	Sodium Thiosulfate Solution (5%)	125 ml	18-25°C
SSC-NFS125	Nuclear Fast Red Solution	125 ml	18-25°C

Precautions: Avoid contact with skin and eyes.
 Harmful if swallowed.
 Follow all Federal, State, and local regulations regarding disposal.
 Use in chemical fume hood whenever possible.

Mixed Storage Conditions. Separate Contents.

For information regarding ordering individual components, please contact us at: 800-442-3573.
Control Slides Available. Catalog: CS-ARG/25, Argentaffin, 25/pack

Preparation of Reagent Prior to Beginning:

Prepare **Ammoniacal Silver Solution** immediately prior to use.

In new or chemically cleaned glassware, mix 27ml Distilled/Deionized water with one vial of Silver Nitrate Solution (10%) and blend completely. Carefully add Concentrated Ammonium Hydroxide (Not included) one drop at a time, swirling gently after each drop. Initially the mixture will turn dark brown and then gradually become transparent with a fine layer of sediment. The solution is ready for immediate use when all sediment dissolves.

Procedure (Standard):

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Place freshly mixed Ammoniacal Silver Solution in a 58-60°C waterbath and allow adequate time for temperature to equalize.
3. Incubate slide in warmed Ammoniacal Silver Solution for 30-60 minutes or until tissue section becomes yellowish/brown in color. (**NOTE:** Melanin typically stains in 30 minutes while Argentaffin stains in 50-60 minutes)
4. Rinse in 3 changes of distilled water.
5. Incubate slide in Gold Chloride Solution (0.2%) for 30 seconds.
6. Rinse in 3 changes of distilled water.
7. Incubate slide in Sodium Thiosulfate Solution (5%) for 1-2 minutes.
8. Rinse for 2 minutes in running tap water followed by 2 changes of distilled water.
9. Incubate slide in Nuclear Fast Red Solution for 5 minutes.
10. Rinse for 2 minutes in running tap water followed by 2 changes of distilled water.
11. Dehydrate very quickly in 3 changes of absolute alcohol.
12. Clear, and mount in synthetic resin.

References:

1. Sheenan, D.C., Hrapchak, B.B. Theory and Practice of Histotechnology, 2nd Edition. Battelle Press, Columbus, OH.
2. Gaitanis, G., et al. Novel application of the Masson-Fontana Stain for Demonstrating Malassezia Species Melanin-Like Pigment Production In Vitro and in Clinical Specimens. Journal of Clinical Microbiology. 2005, August; 43(8): pages 4147-4151.
3. Kimura, M., et al. Fontana-Masson – stained tissue from culture-proven mycoses. Archives of Pathology & Laboratory Medicine. 1998, December; 122(12): page 11.
4. Lazcano, O., et al. Combined Fontana-Masson-Mucin staining of Cryptococcus neoformans. Archives of Pathology & Laboratory Medicine. 1991, November; 115(11): pages 1145-1149.
5. Ro, J.Y., et al. Advantage of Fontana-Masson stain in capsule-deficient cryptococcal infection. Archives of Pathology & Laboratory Medicine. 1987, January; 111(1): pages 53-57.



Instructions For Use **IFU-047** **SSK-FONTMASS**

Rev. Date: Aug. 17, 2016 **Revision: 34** Page 3 of 3

Lot-to-Lot Validation Form Fontana-Masson Stain Kit Catalog: SSK-FONTMASS

Kit Lot Number: _____
Kit Expiration Date: _____
Date Tested: _____
Control Tissue (#) _____
Approved for Use: Y/N _____
Date put into use: _____
If not approved,
corrective actions
taken: _____
Approved by: _____

Kit Component	Lot #
Gold Chloride (0.2%) Sol.	_____
Silver Nitrate (10%) Sol.	_____
Sodium Thiosulfate (5%)	_____
Nuclear Fast Red Solution	_____

Replacement Component if used	Replacement Date	Lot #	Accepted Y/N	Comments
Gold Chloride (0.2%)				
Silver Nitrate Sol. (10%)				
Sodium Thiosulfate (5%)				
Nuclear Fast Red Sol.				
Approved By:				

StatLab is providing this form to assist with reagent lot validation as stated in CLIA'88 Standard 493.1256-For reagent(s), the laboratory must do the following: Check each batch (prepared in-house), lot number (commercially prepared) and shipment of reagents, stains, and identification systems (systems using two or more substrates or two or more reagents, or a combination) when prepared or opened for positive and negative reactivity, if applicable.