

Ultra High Def™ Blue-HRP

CRM-003-30
CRM-003-110

Document #: IFU-CRM-003-Ultra High Def™ Blue-HRP

Release Date: 12/04/2017, IFU-159 Rev A

Intended Use

For In Vitro Diagnostic Use

Summary and Explanation

Ultra High Def™ Blue-HRP is a substrate-chromogen system designed to be used for either IHC or ISH when utilizing horseradish peroxidase. Ultra High Def™ Red-HRP can be permanently mounted to produce a strong azure blue color that can be easily distinguished from other stains.

Principles of the Procedures

Substrate/chromogen in conjunction with horseradish peroxidase (HRP)-based immunostaining or in situ hybridization systems. With working stability of up to 2 hours, this substrate/chromogen system can be used in any automation system as well as manual use.

Reagents Provided

Kit Contents	30 mL	110 mL
UHD™ Blue HRP Substrate Buffer	30 mL	110 mL
UHD™ Blue HRP Chromogen	1 mL	3 mL
Empty Mixing Bottle	1	1

Prepare the Following Solutions Before Use

1. Aliquot 1mL of Ultra High Def™ Blue HRP Substrate Buffer in a mixing bottle.
2. Add one drop (~20µl) of concentrated Ultra High Def™ Blue HRP Chromogen solution.
3. Replace tip, mix, and allow solution to reach room temperature before using.
4. The Ultra High Def™ Red HRP chromogen-substrate working solution is light sensitive and should be kept away from light as much as possible.
5. Working solution is stable for up to 6 hours, for optimal results prepare fresh reagent.

Materials Required But Not Provided

All the reagents and materials required for IHC are not provided. Pretreatment reagents, detection systems, control slides, control reagents and other ancillary reagents are available from StatLab. Please refer to our

website at: www.statlab.com

Storage and Handling

Store at 2-8°C. Do not use product after the expiration date printed on vial. If reagents are stored under a condition other than those specified in the package insert, they must be verified by the user. Diluted reagents should be used promptly.

Staining Procedure

1. Once sections have been incubated with peroxidase, wash sections with wash buffer.
 - a. **Pre-Mix Working Solution: (Automation)** Ultra High Def™ Red HRP has a working solution stability of up to 6 hours and can be loaded directly onto instrument as a single solution. Reduce exposure to light to achieve optimal staining. Working solution is applied directly to slide. Incubate for 10 - 20 minutes.
 - b. **On Board Mixing: (Automation)** Instruments that have on-board mixing capability can load the chromogen and substrate-buffer components independently. Working solution is made mixing reagents 1:50 in on-board mixing station before application to slide. Incubate for 10 - 20 minutes.
 - c. **Manual Use:** Mix substrate-chromogen and buffer in a 1:50 ratio and apply directly to slide. Incubate for 10 – 20 minutes.
2. Counterstain with Hematoxylin or other counterstain. Wash with DI water followed by Immuno wash buffer.
3. Slides should be air dried (do not dehydrate in alcohol). After rinsing off counterstain in DI H₂O. Drain off fluid without further rinsing. Leave slides on benchtop for at least 20 minutes to air dry, and then permanently mount.

Precautions

1. Consult local and/or state authorities with regard to recommended method of disposal.
2. Materials of human or animal origin should be handled as biohazardous materials and disposed of with proper precautions.
3. Avoid microbial contamination of reagents. Contamination could produce erroneous results.
4. This reagent may cause irritation. Avoid contact with eyes and mucous membranes.
5. If reagent contacts these areas, rinse with copious amounts of water.
6. Do not ingest or inhale any reagents

Troubleshooting

If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem is suspected, contact StatLab Medical Supply Headquarters: 2090 Commerce Drive, McKinney, TX 75069. Email our team at ihctech@statlab.com/Call at (800) 442-3573.

