


StatLab Quantum HDx Tubing Cleaning Kit

- Contents:** Tubing Cleaning Component #1-Qty-1 Bottle
 Tubing Cleaning Component #2-Qty-1 Bottle
- Description:** StatLab Quantum HD Tubing Cleaning kit removes diaminobenzidine (DAB) stains and precipitate from the reagents aspiration tubing (Z1 tubing) and helps remove wax residues in the waste tubing (Z2 tubing) of the StatLab Quantum HDx™ IHC Stainer.
- Precautions:** These solutions are harmful and irritating to the eyes, respiratory system and skin. It may cause lung and stomach damage if ingested. Wear appropriate Personal Protective Equipment, including: acid resistant gloves, safety glasses, and laboratory coat or clothing cover. SDS is available.
 Unused solution should be disposed of according to local, State and Federal regulations.
- There are no express or implied warranties which extend beyond this datasheet. StatLab is not liable for personal injury, property damage or economic loss caused by this product.
- Storage:** Store at room temperature. Do not use after date of expiration.
- Instructions for Use:** It is recommended that the tubing cleaning is performed on a routine basis.
1. Load the reagent rack onto the instrument.
 2. Remove the caps and place the Tubing Cleaning Component #1 on reagent rack, position #1 and place Tubing Cleaning Component #2 on reagent rack at position #5. Ensure any caps or lids are removed.



Instructions for Use continued

3. Connect the Wash Buffer in-line. Be sure there is sufficient buffer in the bottle.
4. Log on to the main user interface. Go to System Utilities>Tools and click "Tubing Clean" function.



5. The Z1 and Z2 probes will both aspirate 5 mL of cleaning solution into the tubing from their respective reagent rack positions. The system will then initiate a 20 minute countdown.
6. After 20 minutes has elapsed, the system will automatically perform a system initialization and purge the cleaning solution waste. *The user may click "Cancel" at any moment to immediately purge the waste from the tubing.*
7. Repeat the Tube Cleaning function as needed, to insure there are no visible wax residues and that DAB precipitate or build-up has been removed.

