Dako Gill's 3 Hematoxylin

Code CS708 / Code CS708 / N.º de catálogo CS708

ENGLISH

Intended use	For in vitro diagnostic use. Dako Gill's 3 Hematoxylin is intended for automated primary hematoxylin & eosin (H&E) staining together with Dako Eosin Y Phloxine B (Code CS710), Dako Bluing Buffer (Code CS702) and Dako Differentiation Solution (Code CS707) on the Dako CoverStainer instrument. This ready-to-use hematoxylin is intended for laboratory use as a qualitative histologic stain to demonstrate nuclei (blue) in formalin-fixed, paraffin-embedded (FFPE) tissue sections, frozen tissue sections, and cell preparations. The histochemical stain should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Hematoxylin is extracted from the heartwood of the logwood tree, <i>Haematoxylon campechianum</i> , a species native to Central America (1). The oxidation product hematein forms a cationic (positively charged) complex with a metallic salt (mordant) which enables chelation with negatively charged binding sites of cell components, which are most abundant on chromatin located in the nucleus. Bluing converts the initial red color of the hematoxylin to a blue color. As a result, cell nuclei are stained a blue/purple color (2-4). Refer to our <i>General Instructions for Primary Staining</i> (PD00852) for: Principle of Procedure; Materials Required, but Not Supplied; Storage; Specimen Preparation; Staining Procedure; Quality Control; Staining Evaluation; General Limitations; Troubleshooting. Available at agilent.com/library/eifu .
Reagent provided	Aqueous solution of hematoxylin. The bottle contains 1000 mL ready-to-use reagent. When Dako Gill's 3 is loaded onto the Dako CoverStainer instrument, it must be used within 5 days or replaced after 3000 slides (whichever comes first). The 5-day onboard stability applies when reagent is drained back into the bottle during the night.
Materials required, but not supplied	Dako Eosin Y Phloxine B (Code CS710) Dako Bluing Buffer (Code CS702) Dako Differentiation Solution (Code CS707) Refer to our <i>General Instructions for Primary Staining</i> (PD00852).
Precautions	1. For in vitro diagnostic use. 2. For professional users only. 3. For prescription use only (Rx Only). 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Reagents should be disposed of in accordance with all local, regional, national and international regulations. 6. Safety data sheet is available on www.agilent.com or on request. 7. Dako Gill's 3 Hematoxylin is labeled:

Warning

H226	Flammable liquid and vapour.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
P280	Wear protective gloves. Wear protective clothing. Wear eye or face protection.
P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P241	Use explosion-proof electrical, ventilating or lighting equipment.
P362 + P364	Take off contaminated clothing and wash it before reuse.
P403 + P235	Store in a well-ventilated place. Keep cool.
P501	Dispose of contents and container in accordance with all local, regional, national and international regulations.

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8. Contact Agilent Pathology Support via www.agilent.com to report any unusual results. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the country in which the user and/or the patient is established.

Storage

Store at 15–30 °C. Store in the dark. Do not use after expiration date stamped on bottle. When Dako Gill's 3 is loaded onto the Dako CoverStainer instrument, it must be used within 5 days or replaced after 3000 slides (whichever comes first). The 5-day onboard stability applies when reagent is drained back into the bottle during the night. If reagents are stored under any conditions other than those specified, conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, a daily routine quality control slide(s) containing suitable tissue fixed and processed in a similar manner to the test specimens should be stained prior to routine use to ensure Dako Gill's 3 Hematoxylin is performing as intended. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact Agilent Pathology Support.

Specimen preparation

Paraffin sections:
Dako Gill's 3 Hematoxylin can be used for primary staining of paraffin-embedded tissue sections fixed in formalin.

Frozen sections and cell preparations:
Dako Gill's 3 Hematoxylin can be used for primary staining of acetone-fixed, frozen sections or fixed cell preparations.

Staining procedure

Dako Gill's 3 Hematoxylin (Code CS708) is ready-to-use. For automated use on the Dako CoverStainer instrument. When Dako Gill's 3 is loaded onto the Dako CoverStainer instrument, it must be used within 5 days or replaced after 3000 slides (whichever comes first). The 5-day onboard stability applies when reagent is drained back into the bottle during the night. Please read the User Guide for the Dako CoverStainer instrument and program the instrument accordingly with the regressive H&E staining protocol provided below.

Dako Regressive H&E protocol**:

Pre-heating for 10 minutes at 65 °C

Step	Dip tank	Station	Bottle	Reagent	Time min	Time max
1	1	1-2	A	Xylene 1*	03:00	03:30
2	2	3-4	B	Xylene 2*	03:00	03:30
3	3	5-6	C	96 % EtOH 1	01:00	-
4	4	7-8	D	96 % EtOH 2	01:00	-
5	5	9-10	E	70 % EtOH 1	02:00	-
6	18	35-36	-	Tap water	01:00	-
7	8	15	H	Dako Gill's 3 Hematoxylin	01:45	01:45
8	7	13-14	G	Deionized water 1	01:00	01:00
9	9	18	K	Differentiation Solution	00:30	00:30
10	17	33-34	-	Tap water	01:00	01:00
11	9	17	J	Dako Bluing Buffer	01:00	01:00
12	17	33-34	-	Tap water	01:00	01:00
13	6	11-12	F	70 % EtOH 2	01:00	01:00
14	10	19	L	Dako Eosin Y Phloxine B	02:00	02:00
15	11	21	N	96 % EtOH 3	01:00	01:00
16	11	22	O	99.9 % EtOH 1	01:00	01:00
17	12	23-24	P	99.9 % EtOH 2	01:00	01:00
18	13	25-26	Q	99.9 % EtOH 3	01:00	01:00
19	14-16	27-32	R, S, T	Exit solvent*	01:00	-

*Xylene or xylene substitute, Histo-Clear II (National Diagnostic, Inc), HS-202, is recommended.

**The user must always read the Instructions for Use for the reagents used and consult the Dako CoverStainer User Guide for details. Protocol settings must be verified before use.

After the staining procedure has been completed, the sections are mounted onboard Dako CoverStainer instrument using Dako Mounting Medium (Code CS703) or Dako Toluene-Free Mounting Media (Code CS705). The slides must have a drying period of 10 minutes on the instrument prior to evaluation. Refer to our *General Instructions for Primary Staining* (PD00852).

Quality control

A daily routine quality control slide(s) containing suitable tissue fixed and processed in a similar manner to the test specimens should be stained prior to routine use to ensure Dako Gill's 3 Hematoxylin is performing as intended.

Staining pattern

Dako Gill's 3 Hematoxylin stains cell nuclei blue.

Additional information

The user should read the User Guide for the Dako CoverStainer instrument and product-specific Instructions for Use carefully and become familiar with all the components and the instrumentation prior to use.

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Agilent Technologies Singapore (International) Pte Ltd.
*Xileno o sustituto de xileno. Se recomienda Histo-Clear II (National Diagnostic, Inc), HS-202.