Human Colon Tumor

Colloidal Iron Stain Kit Procedure

100ml Kit Item #: KTCIR | **Liter Kit Item#:** N/A
Pint Kit Item #: KTCIRPT | **Gallon Kit Item#:** N/A

Control Slide(s)	Item#	Included Components	
Colon	CSC1025P	Colloidal Iron Stock	12% Acetic Acid
Small Intestine	CSS0325P	3% Hydrochloric Acid	Van Gieson's Stain
		3% Potassium Ferrocyanide	
Components Not Included:			
Glacial Acetic Acid (Item# SPA0243)			

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify acid mucopolysaccharides. Acid mucopolysaccharides stain bright blue and collagen shades of red.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.


WORKING COLLOIDAL IRON SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	15ml	Distilled Water	Into a coplin jar or container.
2	Add	5ml	Glacial Acetic Acid (<i>Not Included</i>)	Mix thoroughly.
3	Add	20ml	Colloidal Iron Stock	Mix thoroughly.

WORKING IRON STAIN SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	20ml	3% Potassium Ferrocyanide	Into a coplin jar or container.
2	Add	20ml	3% Hydrochloric Acid	Mix thoroughly.

STAINING PROCEDURE:

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	12% Acetic Acid	--	--	30	Once complete, discard solution.
5	Immerse	Working Colloidal Iron Solution	--	30	--	Agitate periodically.
6	Rinse	12% Acetic Acid, 3 changes	--	2	--	2 minutes each change. Once complete, discard solution.
7	Immerse	Working Iron Stain Solution	--	10	--	Agitate periodically. Once complete, rinse in running DI water (1 minute) and continue.
8	Immerse	Van Gieson's Stain	--	--	30-45	Without rinsing, continue to next step.
9	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

- Muller G: ACTA Histochem (Jena); 2:68, 1955.
- With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Working Colloidal Iron Solution and Standard Working Iron Stain Solution preparation procedures yields approximately 40ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 8 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 19 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	100	900	9	250ml Glass Stain Dish	8	675	85
30ml Glass Coplin Jar	66	891	14	200ml Bath Autostainer	10	720	72
40ml Hellendahl Jar	50	900	18	400ml Bath Autostainer	5	765	153

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS07	Irritant	CE

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Muller G: ACTA Histochem (Jena); 2:68, 1955.
2. With modifications by AMTS R&D Department, 1979-2019.

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