Human Aorta

Verhoeff's Elastic Stain Kit Procedure

100ml Kit Item #: KTVEL | **Liter Kit Item#:** KTVELLT
Pint Kit Item #: KTVELPT | **Gallon Kit Item#:** N/A

Control Slide(s)	Item#	Included Components	
Aorta	CSA0825P	5% Alcoholic Hematoxylin	5% Sodium Thiosulfate
		10% Ferric Chloride	Van Gieson's Stain
		Universal Iodine Solution™	

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify elastic fibers. Elastic fibers stain black, collagen red, and other tissue yellow.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-6µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label. Some permanent Mounting Mediums may cause fading over time with picric acid based stains. While MicroMount™ or PermaSlip™ is recommended for optimal staining preservation, many other mounting mediums may also provide successful results. Be sure to test your preferred Mounting Medium with this stain if not using one of the recommended products.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.


VERHOEFF'S ELASTIC STAIN PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	40ml	5% Alcoholic Hematoxylin	Into a chemically cleaned container or new/unused plasticware.
2	Add	20ml	10% Ferric Chloride	Mix thoroughly
3	Add	20ml	Universal Iodine Solution™	Mix thoroughly.
4	Filter	--	Paper Towel	Filter into a chemically cleaned container or new/unused plasticware.

2% FERRIC CHLORIDE PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	40ml	Distilled/DI Water	Into a chemically cleaned container or new/unused plasticware.
2	Add	10ml	10% Ferric Chloride	Mix thoroughly

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Verhoeff's Elastic Stain	--	15	--	Once complete, rinse in running tap water (1 minute) and continue.
5	Immerse	2% Ferric Chloride	--	1-5	--	Differentiate until elastic fibers are black (use microscope). Once complete, rinse in running tap water (1 minute) and continue.
6	Immerse	5% Sodium Thiosulfate	--	1	--	Once complete, rinse in running tap water (1-2 minutes) and continue.
7	Immerse	Van Gieson's Stain	--	--	30	
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	See storage and usage notes below.

***NOTE:** Some permanent Mounting Mediums may cause fading over time with picric acid based stains. While MicroMount™ or PermaSlip™ is recommended for optimal staining preservation, many other mounting mediums may also provide successful results. Be sure to test your preferred Mounting Medium with this stain if not using one of the recommended products.

1. Mallory: Pathological Techniques, pg 170.
2. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Verhoeff's Working Elastic Stain Solution preparation procedure yields approximately 80ml of solution and Standard 2% Ferric Chloride Solution preparation procedure yields approximately 50ml of solution and both must be scaled up to accommodate desired autostainer bath size. A minimum of 5 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 16 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	225	9	250ml Glass Stain Dish	2	169	18
30ml Glass Coplin Jar	16	216	14	200ml Bath Autostainer	2	144	72
40ml Hellendahl Jar	12	216	18	400ml Bath Autostainer	1	153	153

CE MARKINGS AND DESIGNATIONS:

	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
	Batch Code		Use By		Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Flammable	
	Corrosive		Irritant		Health Hazard	

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Mallory: Pathological Techniques, pg 170.
2. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 5% Alcoholic Hematoxylin, 10% Ferric Chloride, Universal Iodine Solution™, 5% Sodium Thiosulfate, Van Gieson's Stain

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