Normal Human Colon

## Alcian Blue pH 2.5 Stain Kit Procedure

100ml Kit Item #: KTABP2.5 | Liter Kit Item#: KTABP2.5LT  
 Pint Kit Item #: KTABP2.5PT | Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components
Colon	CSC1025P	3% Acetic Acid
Small Intestine	CSS0325P	Alcian Blue Stain pH 2.5 Nuclear Fast Red Stain


**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify acid mucins, carboxylated and sulfated mucosubstances. Strongly sulfated mucosubstances are not demonstrated. Acid mucins, carboxylated and sulfated mucosubstances are stained blue, cytoplasm and background pale pink, and nuclei pink to red.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

### STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	3% Acetic Acid	--	3	--	Without rinsing, continue to next step.
5	Immerse	Alcian Blue Stain pH 2.5	--	30	--	Once complete, rinse in running tap water (1 minute) and continue.
6	Immerse	Nuclear Fast Red Stain	--	5	--	Once complete, rinse in running tap water (1 minute) and continue.
7	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
8	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
9	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology, 1980, 172 – 173
2. With modifications by AMTS R&D Department, 1979-2019.

**AUTOSTAINER CONFIGURATION AND NOTES:**

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 4 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 15 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	500	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	16	480	30	200ml Bath Autostainer	2	320	160
40ml Hellendahl Jar	12	480	40	400ml Bath Autostainer	1	340	340

**CE MARKINGS AND DESIGNATIONS:**

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By	<b>EC REP</b>	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Health Hazard	<b>CE</b>

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

**StatLab**  
2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology, 1980, 172 – 173
2. With modifications by AMTS R&D Department, 1979-2019.

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**COMPOSANTS INCLUS:** 3% Acetic Acid, Alcian Blue Stain pH 2.5, Nuclear Fast Red Stain

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