Human Brain

## Modified Bielschowsky's Stain Kit Procedure

**100ml Kit Item #:** KTBIE

**Liter Kit Item#:** N/A

**Pint Kit Item #:** N/A

**Gallon Kit Item#:** N/A

| Control Slide(s) | Item#    |
|------------------|----------|
| Alzheimer        | CSA0225P |

| Included Components        |                            |
|----------------------------|----------------------------|
| Silver Nitrate Solution    | Stock Developer Solution C |
| Stock Developer Solution A | 5% Sodium Thiosulfate      |
| Stock Developer Solution B |                            |

**Components Not Included:**  
Concentrated Ammonium Hydroxide (Item#: SPA1326)

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify axons, intracellular neurofibrils and senile plaques. Intracellular neurofibrils (neurofibrillary tangles), axons, and senile plaques stain gray to black, and nuclei and background yellow to gold.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 6-8µm tissue sections or frozen sections.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

**AMMONIACAL SILVER SOLUTION PREPARATION:** Prepare solution at time of use. Solution expires after one use.

| # | Action | Amount | Chemical/Reagent                               | Details  |
|---|--------|--------|--|--|
| 1 | Add    | 20ml   | Silver Nitrate Solution                        | Into a chemically clean container or new/unused plasticware.   |
| 2 | Add    | Drops  | Concentrated Ammonium Hydroxide (Not Included) | Add drop by drop while carefully agitating container until <u>no precipitate remains</u> . Solution will turn brown, then colorless. |

**DEVELOPER SOLUTION PREPARATION:** Prepare solution at time of use. Solution expires after one use.

| # | Action | Amount  | Chemical/Reagent             | Details  |
|---|--------|---------|------------------------------|--|
| 1 | Add    | 50ml    | Distilled/DI Water           | Into a chemically clean container or new/unused plasticware. |
| 2 | Add    | 8 Drops | Stock Developer Solution "A" | Mix thoroughly.  |
| 3 | Add    | 8 Drops | Stock Developer Solution "B" | Mix thoroughly.  |
| 4 | Add    | 4 Drops | Stock Developer Solution "C" | Mix thoroughly.  |

**AMMONIA WATER PREPERATION:** Prepare solution at time of use.

| # | Action | Amount  | Chemical/Reagent                               | Details  |
|---|--------|---------|--|--|
| 1 | Add    | 25ml    | Distilled/DI Water                             | Into a chemically clean container or new/unused plasticware. |
| 2 | Add    | 3 Drops | Concentrated Ammonium Hydroxide (Not Included) | Mix thoroughly.  |

**STAINING PROCEDURE:** Preheating Required. See steps 4 and 5 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

| #  | Action        | With                            | Heat °C | Time |      | Details   |
|----|---------------|---------------------------------|---------|------|------|---|
|    |               |                                 |         | Mins | Secs |   |
| 1  | Deparaffinize | Xylene or Substitute, 2 changes | --      | 5    | --   | 5 minutes each change or as required if using a xylene substitute.  |
| 2  | Rinse         | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change or as required if using graded alcohols.   |
| 3  | Rinse         | Running DI Water                | --      | 2    | --   |   |
| 4  | Waterbath     | Silver Nitrate Solution         | 40°     | 15   | --   | Immerse into preheated solution. Once complete, <b>rinse in running DI water (2 minutes)</b> and continue.            |
| 5  | Waterbath     | Ammoniacal Silver Solution      | 40°     | 10   | --   | Immerse into preheated solution. Once complete, drain excess solution. Without rinsing, continue to next step         |
| 6  | Immerse       | Developer Solution              | --      | --   | 3-30 | Agitate. Remove when golden brown tissue color is achieved.   |
| 7  | Immerse       | Ammonia Water                   | --      | --   | 30   | Quickly immerse to stop developing process. Once complete, <b>rinse in running DI water (2 minutes)</b> and continue. |
| 8  | Immerse       | 5% Sodium Thiosulfate           | --      | 2    | --   | Once complete, <b>rinse in running DI water (2 minutes)</b> and continue.   |
| 9  | Dehydrate     | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change.   |
| 10 | Clear         | Xylene or Substitute, 3 changes | --      | 1    | --   | 1 minutes each change or as required if using a xylene substitute.  |
| 11 | Coverslip     | Permanent Mounting Media        | --      | --   | --   |   |

1. Mirra SS Hart MN Terry RD: Making the Diagnosis of Alzheimer's Disease, Arch Pathol Lab Med. 1993:117: 132 - 144.
2. Luna LG: Histopathologic Methods and Color Atlas of Special Stains and Tissue Artifacts; American Histolabs, Inc. 1992: 483 - 485.
3. With modifications by AMTS R&D Department, 1979-2019.

**AUTOSTAINER CONFIGURATION AND NOTES:**

This stain kit is not available in the pint or larger size.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

| Bath Type              | Uses | Slides | S.C. | Bath Type              | Uses | Slides | S.C. |
|------------------------|------|--------|------|------------------------|------|--------|------|
| 20ml Plastic Slide Jar | 12   | 48     | 4    | 250ml Glass Stain Dish | n/a  | n/a    | n/a  |
| 30ml Glass Coplin Jar  | 8    | 64     | 8    | 200ml Bath Autostainer | n/a  | n/a    | n/a  |
| 40ml Hellendahl Jar    | 6    | 96     | 16   | 400ml Bath Autostainer | n/a  | n/a    | n/a  |

**CE MARKINGS AND DESIGNATIONS:**

|            |                                    |           |                                   |               |                          |  |
|------------|------------------------------------|-----------|-----------------------------------|---------------|--------------------------|--|
| <b>REF</b> | Catalogue Number                   |           | Temperature Limitation            |               | Manufacturer             | American MasterTech Scientific<br>1330 Thurman St.<br>Lodi, CA 95240 USA<br>Tel 800 860 4073<br>Fax 209 368 4136 |
| <b>LOT</b> | Batch Code                         |           | Use By                            | <b>EC REP</b> | Representative           | Emergo Europe<br>Prinsessegracht 20<br>2514 AP The Hague<br>The Netherlands                                      |
| <b>IVD</b> | In Vitro Diagnostic Medical Device |           | Consult Instructions Prior to Use | <br>GHS05     | Corrosive                | <b>CE</b>  |
| <br>GHS07  | Irritant                           | <br>GHS08 | Health Hazard                     | <br>GHS09     | Environmentally Damaging |  |

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

**StatLab**  
2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Mirra SS Hart MN Terry RD: Making the Diagnosis of Alzheimer's Disease, Arch Pathol Lab Med. 1993;117: 132 - 144.
2. Luna LG: Histopathologic Methods and Color Atlas of Special Stains and Tissue Artifacts; American Histolabs, Inc. 1992: 483 - 485.
3. With modifications by AMTS R&D Department, 1979-2019.

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