Animal Liver

Chandler's Precision® Reticulum Stain Kit Procedure

100ml Kit Item #: KTCPR

Liter Kit Item#: KTCPRLT

Pint Kit Item #: KTCPRPT

Gallon Kit Item#: N/A

| Control Slide(s) | Item# |
|------------------|----------|
| Reticulum | CSR0125P |

| Included Components | |
|----------------------------|------------------------|
| Type 1 Silver Solution | 20% Formalin |
| 3% Sodium Hydroxide | 0.1% Gold Chloride |
| 1% Potassium Permanganate | 5% Sodium Thiosulfate |
| 3% Potassium Metabisulfite | Nuclear Fast Red Stain |
| 3% Ferric Ammonium Sulfate | |

Components Not Included:

Concentrated Ammonium Hydroxide (Item#: SPA1326)

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify reticulum fibers. Reticulum fibers stain gray to black and nuclei pink to red.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

AMMONIACAL SILVER SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use. Ammoniacal Silver Solution reacts with metals. Do not use metal forceps, racks, or tools with the solution.

| # | Action | Amount | Chemical/Reagent | Details |
|---|--------|--------|---|---|
| 1 | Add | 1 vial | Type 1 Silver Solution | Into a chemically cleaned container or new/unused plasticware. |
| 2 | Add | 17ml | Distilled Water | Mix thoroughly. |
| 3 | Add | 1ml | Concentrated Ammonium Hydroxide (Not Included) | Drop by Drop. Mix thoroughly; solution will turn brown, then colorless. |
| 4 | Add | 10ml | 3% Sodium Hydroxide | Mix thoroughly. If no precipitate remains, skip to Step 6. |
| 5 | Add | Drops | Concentrated Ammonium Hydroxide (Not Included) | Drop by Drop while mixing, until no precipitate remains. |
| 6 | Add | 30ml | Distilled Water | Mix thoroughly. |

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

| # | Action | With | Heat °C | Time | | Details |
|----|---------------|---------------------------------|---------|------|-------|--|
| | | | | Mins | Secs | |
| 1 | Deparaffinize | Xylene or Substitute, 2 changes | -- | 5 | -- | 5 minutes each change or as required if using a xylene substitute. |
| 2 | Rinse | Absolute Alcohol, 3 changes | -- | 1 | -- | 1 minute each change or as required if using graded alcohols. |
| 3 | Rinse | Running DI Water | -- | 1 | -- | |
| 4 | Immerse | 1% Potassium Permanganate | -- | 5-10 | -- | Once complete, rinse in running DI water (1 minute) and continue. |
| 5 | Immerse | 3% Potassium Metabisulfite | -- | -- | 30-60 | Until section is colorless. Once complete, rinse in running DI water (1 minute) and continue. |
| 6 | Immerse | 3% Ferric Ammonium Sulfate | -- | 10 | -- | Once complete, rinse in running DI water (1 minute) and continue. |
| 7 | Immerse | Ammoniacal Silver Solution | -- | 2-3 | -- | Once complete, rinse in running DI water (1 minute) and continue. |
| 8 | Immerse | 20% Formalin | -- | 1 | -- | Once complete, rinse in running DI water (1 minute) and continue. |
| 9 | Immerse | 0.1% Gold Chloride | -- | 3-5 | -- | Once complete, rinse in running DI water (1 minute) and continue. |
| 10 | Immerse | 5% Sodium Thiosulfate | -- | 1-2 | -- | Once complete, rinse in running DI water (1 minute) and continue. |
| 11 | Immerse | Nuclear Fast Red Stain | -- | 3-5 | -- | Once complete, rinse in running DI water (1 minute) and continue. |
| 12 | Dehydrate | Absolute Alcohol, 3 changes | -- | 1 | -- | 1 minute each change. |
| 13 | Clear | Xylene or Substitute, 3 changes | -- | 1 | -- | 1 minute each change or as required if using a xylene substitute. |
| 14 | Coverslip | Permanent Mounting Media | -- | -- | -- | |

- Gomori: A Modification of the Silver Impregnation Method of Staining Reticular Fibers; American Journal of Clinical Pathology, 21:1951, 897-899
- With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Ammoniacal Silver Solution preparation procedure yields approximately 61ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 9 baths are required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 20 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

| Bath Type | Uses | Slides | S.C. | Bath Type | Uses | Slides | S.C. |
|------------------------|------|--------|------|------------------------|------|--------|------|
| 20ml Plastic Slide Jar | 60 | 240 | 4 | 250ml Glass Stain Dish | 4 | 120 | 30 |
| 30ml Glass Coplin Jar | 40 | 320 | 8 | 200ml Bath Autostainer | 6 | 180 | 30 |
| 40ml Hellendahl Jar | 30 | 480 | 16 | 400ml Bath Autostainer | 3 | 150 | 50 |

CE MARKINGS AND DESIGNATIONS:

| | | | | | | |
|------------|------------------------------------|-----------|-----------------------------------|---------------|--------------------------|--|
| REF | Catalogue Number | | Temperature Limitation | | Manufacturer | American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136 |
| LOT | Batch Code | | Use By | EC REP | Representative | Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands |
| IVD | In Vitro Diagnostic Medical Device | | Consult Instructions Prior to Use | GHS05 | Corrosive | CE |
| GHS07 | Irritant | GHS08 | Health Hazard | GHS09 | Environmentally Damaging | |

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

- Gomori: A Modification of the Silver Impregnation Method of Staining Reticular Fibers; American Journal of Clinical Pathology, 21:1951, 897-899
- With modifications by AMTS R&D Department, 1979-2019.

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