Human Adrenal Gland

Congo Red Amyloid Stain Kit Procedure

100ml Kit Item #: KTCRE

Liter Kit Item#: KTCRELT

Pint Kit Item #: KTCREPT

Gallon Kit Item#: N/A

Control Slide(s)

Amyloid Alpha

Item#

CSA0425P

Included Components

Congo Red Amyloid Stain

1% Sodium Hydroxide

Modified Mayer's Hematoxylin

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens or frozen sections (in vitro) to identify amyloid deposits. Amyloid stain pale pink to red (Birefringent if viewed using a polarized light microscope) and nuclei blue.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded or frozen tissue section cut at 8-10µm.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Congo Red Amyloid Stain	--	60	--	Without rinsing, continue to next step.
5	Immerse	1% Sodium Hydroxide	--	--	15-20	Differentiate until Amyloid is pink to red. Use a quick tap water rinse and microscope to check for correct differentiation. Once complete, rinse in running tap water (1 minute) and continue.
6	Immerse	Modified Mayer's Hematoxylin	--	3	--	Once complete, rinse in running tap water (1 minute) and continue.
7	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
8	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
9	Coverslip	Permanent Mounting Media	--	--	--	

1. A.F.I.P. Manual Laboratory Methods in Histotechnology; pg 168, 1992.
2. Lilly; Histopathologic Technique; pg 293.
3. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 4 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 15 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	50	1000	20	250ml Glass Stain Dish	4	750	188
30ml Glass Coplin Jar	33	990	30	200ml Bath Autostainer	5	800	160
40ml Hellendahl Jar	25	1000	40	400ml Bath Autostainer	2	680	340

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Irritant		Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. A.F.I.P. Manual Laboratory Methods in Histotechnology; pg 168, 1992.
2. Lilly; Histopathologic Technique; pg 293.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Congo Red Amyloid Stain, 1% Sodium Hydroxide, Modified Mayer's Hematoxylin

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COMPONENTES INCLUIDOS: Congo Red Amyloid Stain, 1% Sodium Hydroxide, Modified Mayer's Hematoxylin

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