Human Skin

Chicago Sky Blue Stain Kit Procedure

30ml Kit Item #: KTCSB

Control Slide(s)	Item#
N/A	

Included Components
Potassium Hydroxide Solution
Chicago Sky Blue Stain

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to rapidly diagnose dermatophytic infections, including Dermatophytosis, Candidial Dermatomyces, and Pityriasis Versicolor (in vitro). Dermatophytes are stained blue and tissue colorless to pink.

SPECIMEN CRITERIA: Skin scraping, hair, or nail clipping collected appropriately and placed onto a microscope slide.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information

STAINING PROCEDURE:

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Add	Potassium Hydroxide Solution	--	--	--	Add one drop to slide.
2	Add	Chicago Sky Blue Stain	--	--	--	Add one drop to slide.
3	Coverslip	Coverslip	--	30	--	Apply coverslip and allow to stand at room temperature for 30 minutes. Once complete, review with a microscope.

- Fonseka S, Lim CSH, Bandara U, Dissanayake M: New Contrast Stain for the Rapid Diagnosis of Dermatophytosis and Pityriasis Versicolor. ASCP Lab Medicine, 2011; 42: 649-652
- With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES: N/A

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	n/a	n/a	n/a	250ml Glass Stain Dish	n/a	n/a	n/a
30ml Glass Coplin Jar	n/a	n/a	n/a	200ml Bath Autostainer	n/a	n/a	n/a
40ml Hellendahl Jar	n/a	n/a	n/a	400ml Bath Autostainer	n/a	n/a	n/a

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Corrosive		Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Fonseca S, Lim CSH, Bandara U, Dissanayake M: New Contrast Stain for the Rapid Diagnosis of Dermatophytosis and Pityriasis Versicolor. ASCP Lab Medicine, 2011; 42: 649-652
2. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Potassium Hydroxide Solution, Chicago Sky Blue Stain

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1. Fonseka S, Lim CSH, Bandara U, Dissanayake M: New Contrast Stain for the Rapid Diagnosis of Dermatophytosis and Pityriasis Versicolor. *ASCP Lab Medicine*, 2011; 42: 649-652
2. With modifications by AMTS R&D Department, 1979-2019.

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-
1. Fonseka S, Lim CSH, Bandara U, Dissanayake M: New Contrast Stain for the Rapid Diagnosis of Dermatophytosis and Pityriasis Versicolor. ASCP Lab Medicine, 2011; 42: 649-652
 2. With modifications by AMTS R&D Department, 1979-2019.