Human Lymphoid Tissue

Fite's Stain Kit Procedure (Leprosy & Nocardia Procedure)

100ml Kit Item #: KTFIT

Liter Kit Item#: N/A

Pint Kit Item #: KTFITPT

Gallon Kit Item#: N/A

Control Slide(s) | Item#

Nocardia | CSN0125P

Included Components

1% Acid Alcohol | Advanced Carbol Fuchsin™

Xylene - Peanut Oil | Light Green Counterstain

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify *Mycobacterium leprae* or *Nocardia*. *Mycobacterium leprae* or *Nocardia* stain red and background green.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Filter Advanced Carbol Fuchsin™ when a shiny gold layer appears on the surface. Store/Use each component according to the temperature and expiration on the label. Methylene Blue Stain (AMTS Item# STMBL100, STMBLPT) is available as an alternative counterstain from American MasterTech Scientific.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

DILUTE ACID ALCOHOL: For Nocardia Procedure only! Prepare solution at time of use. Solution expires after one use.

| # | Action | Amount | Chemical/Reagent | Details |
|---|--------|--------------|------------------|--------------------------------|
| 1 | Add | -- | Distilled Water | Into a coplin jar or container |
| 2 | Add | Equal Volume | 1% Acid Alcohol | Mix thoroughly. |

LEPROSY & NOCARDIA STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

| # | Action | With | Heat °C | Time | | Details |
|---|---|--------------------------------|---------|------|------|---|
| | | | | Mins | Secs | |
| 1 | Deparaffinize | Xylene - Peanut Oil, 2 changes | -- | 12 | -- | 12 minutes each change. |
| 2 | Air Dry | On Benchtop, or in Empty Bath | -- | 15 | -- | Remaining oil film will prevent de-staining of <i>Mycobacterium leprae</i> or <i>Nocardia</i> during differentiation. |
| 3 | Rinse | Running Tap Water | -- | 5 | -- | |
| 4 | Immerse | Advanced Carbol Fuchsin™ | -- | 15 | -- | Once complete, rinse in running tap water (1 minute) and continue. |
| 5 | To stain Leprosy follow step 5a. For Nocardia follow step 5b. | | | | | |

| 5a | Leprosy | | | | | |
|---------|-----------------|----|----|-------|----|---|
| Immerse | 1% Acid Alcohol | -- | -- | 30-60 | -- | Differentiate until background is pale pink. Use a quick tap water rinse and microscope to check for correct differentiation. Once complete, rinse in running tap water (1 minute) and continue at step 6. |

| 5b | Nocardia | | | | | |
|---------|---------------------|----|----|---|----|--|
| Immerse | Dilute Acid Alcohol | -- | -- | 1 | -- | 2 dips for 1 second each until background is pale pink. Use a quick tap water rinse and microscope to check for correct differentiation. Once complete, rinse in running tap water (1 minute) and continue at step 6. |

| | | | | | | |
|---|-----------|-------------------------------|----|------|-------|---|
| 6 | Immerse | Light Green Counterstain | -- | -- | 30-60 | To desired intensity. Once complete, rinse in running tap water (10 seconds) and continue. |
| 7 | Air Dry | On Benchtop, or in Empty Bath | -- | 5-10 | -- | Until section is completely dry. |
| 8 | Clear | Xylene or Substitute | -- | -- | 1 | 5-10 dips for 1 second each. |
| 9 | Coverslip | Permanent Mounting Media | -- | -- | -- | |

1. Mallory: Pathological Technique; pg 275.
2. Crowder C Taylor HW: Modified Fite Stain for Demonstration of Mycobacterium Species in Tissue Sections; J of Histotechnology; Vol 19; 2:133-134.
3. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 4 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 8 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

| Bath Type | Uses | Slides | S.C. | Bath Type | Uses | Slides | S.C. |
|------------------------|------|--------|------|------------------------|------|--------|------|
| 20ml Plastic Slide Jar | 25 | 330 | 14 | 250ml Glass Stain Dish | 2 | 248 | 124 |
| 30ml Glass Coplin Jar | 16 | 317 | 20 | 200ml Bath Autostainer | 2 | 211 | 106 |
| 40ml Hellendahl Jar | 12 | 317 | 27 | 400ml Bath Autostainer | 1 | 224 | 225 |

CE MARKINGS AND DESIGNATIONS:

| | | | | | | |
|--|------------------------------------|--|-----------------------------------|---|----------------|--|
| REF | Catalogue Number |  | Temperature Limitation |  | Manufacturer | American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136 |
| LOT | Batch Code |  | Use By | EC REP | Representative | Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands |
| IVD | In Vitro Diagnostic Medical Device |  | Consult Instructions Prior to Use | CE | | |
|  GHS02 | Flammable |  GHS05 | Corrosive |  GHS06 | Toxic | |
|  GHS07 | Irritant |  GHS08 | Health Hazard | | | |

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Mallory: Pathological Technique; pg 275.
2. Crowder C Taylor HW: Modified Fite Stain for Demonstration of Mycobacterium Species in Tissue Sections; J of Histotechnology; Vol 19; 2:133-134.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 1% Acid Alcohol, Xylene Peanut Oil, Advanced Carbol Fuchsin™, Light Green Counterstain

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