Human Skin

Fontana - Masson Stain Kit Procedure

100ml Kit Item#: KTFMA | **Liter Kit Item#:** KTFMALT
Pint Kit Item #: KTFMAPT | **Gallon Kit Item#:** N/A

Control Slide(s) | **Item#**
 Skin | CSS0825P

Included Components

Type 1 Silver Solution	5% Sodium Thiosulfate
0.1% Gold Chloride	Nuclear Fast Red Stain

Components Not Included:
 Concentrated Ammonium Hydroxide (Item#: SPA1326)

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify argentaffin cells and melanin. Melanin and Argentaffin granules stain brown to black and nuclei pink.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

AMMONIACAL SILVER SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use. Ammoniacal Silver Solution reacts with metals. Do not use metal forceps, racks, or tools with the solution.

#	Action	Amount	Chemical/Reagent	Details
1	Add	27ml	Distilled Water	Into a chemically cleaned container or new/unused plasticware.
2	Add	1 vial	Type 1 Silver Solution	Mix thoroughly.
3	Add	Drops	Concentrated Ammonium Hydroxide (Not Included)	Drop by Drop. Mix thoroughly; solution will turn brown, then colorless. Solution is ready to use when all sediments are dissolved.

STAINING PROCEDURE: Preheating Required. See step 4 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Waterbath	Ammoniacal Silver Solution	60°	30-60	--	Immerse into preheated solution. Incubate until section becomes yellowish-brown. Melanin typically stains in approximately 30 minutes; Argentaffin in 50-60 minutes. Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	0.1% Gold Chloride	--	1	--	Until Melanin or Argentaffin granules are brown to black (use microscope). Once complete, rinse in running DI water (1 minute) and continue.
6	Immerse	5% Sodium Thiosulfate	--	1-2	--	Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	Nuclear Fast Red Stain	--	5	--	Once complete, rinse in running DI water (1 minute) and continue.
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 277.
2. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Ammoniacal Silver Solution preparation procedure yields approximately 28ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 5 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 16 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	30	120	4	250ml Glass Stain Dish	2	60	30
30ml Glass Coplin Jar	20	160	8	200ml Bath Autostainer	3	90	30
40ml Hellendahl Jar	15	240	16	400ml Bath Autostainer	1	50	50

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Corrosive	CE
	Irritant		Health Hazard		Dangerous for the Environment	

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Sheehan DC Hrapchak BB; Theory and Practice of Histotechnology; 1980, 277.
2. With modifications by AMTS R&D Department, 1979-2019.

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