Colon

Goldner's Trichrome Stain Kit Procedure

100ml Kit Item #: N/A

Liter Kit Item#: KTGRTL

Pint Kit Item #: N/A

Gallon Kit Item#: N/A

Control Slide(s)

Item#

Trichrome

CST0125P

Included Components

Bouin's Fluid	1% Acetic Acid
Weigert's Hematoxylin A	Phosphomolybdic-Orange G
Weigert's Hematoxylin B	Goldner's Trichrome Light Green Stain
Goldner's Trichrome	
Ponceau-Fuchsin	

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify collagen and muscle. Cytoplasm and muscle are stained red to fuchsia, erythrocytes orange to red, collagen green, and nuclei blue to black.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-6µm tissue section

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

WORKING WEIGERT'S HEMATOXYLIN PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	25ml	Weigert's Hematoxylin A	Into a chemically clean container or new/unused plasticware.
2	Add	25ml	Weigert's Hematoxylin B	Mix thoroughly.

STAINING PROCEDURE: Preheating required. See step 4 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Waterbath	Bouin's Fluid	56-62°	60	--	Immerse into preheated solution. (Optional: Immerse overnight at room temperature). Once complete, rinse in running tap water (3 minutes) until tissue is colorless and continue.
5	Immerse	Working Weigert's Hematoxylin	--	10	--	Once complete, rinse in running tap water (5 minutes) and continue.
6	Immerse	Goldner's Trichrome Ponceau-Fuchsin	--	5	--	
7	Rinse	1% Acetic Acid	--	1	--	
8	Immerse	Phosphomolybdic -Orange G	--	30	--	Until collagen is decolorized.
9	Rinse	1% Acetic Acid	--	1	--	
10	Immerse	Goldner's Trichrome Light Green Stain	--	5-10	--	
11	Rinse	1% Acetic Acid	--	5	--	
12	Dehydrate	95% Reagent Alcohol, 2 changes	--	1	--	1 minute each change.
13	Dehydrate	Absolute Alcohol, 2 changes	--	1	--	1 minute each change.
14	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
15	Coverslip	Permanent Mounting Media	--	--	--	

- Luna, L. G. (1992). Histopathologic methods and color atlas of special stains and tissue artifacts. Gaithersburg, MD: American Histolabs. p 451 – 453
- With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 9 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 21 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	33	297	9	250ml Glass Stain Dish	2	169	85
30ml Glass Coplin Jar	22	297	14	200ml Bath Autostainer	3	216	72
40ml Hellendahl Jar	16	288	18	400ml Bath Autostainer	1	153	153

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS02	Flammable	
 GHS05	Corrosive	 GHS07	Irritant	 GHS08	Health Hazard	

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

- Luna, L. G. (1992). Histopathologic methods and color atlas of special stains and tissue artifacts. Gaithersburg, MD: American Histolabs. p 451 – 453
- With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Bouin's Fluid, Hematoxylin A, Weigert's Hematoxylin B, Goldner's Trichrome Ponceau-Fuchsin, 1% Acetic Acid, Phosphomolybdic-Orange G, Goldner's Trichrome Light Green Stain

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