Human Gall Bladder

## Hall's Bile Stain Kit Procedure

**100ml Kit Item#:** KTHBI

**Liter Kit Item#:** N/A

**Pint Kit Item#:** KTHBIPT

**Gallon Kit Item#:** N/A

Control Slide(s)	Item#
Bile	CSB0125P

Included Components
Fouchet's Reagent
Van Gieson's Stain

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify bile or bilirubin. Bile or Bilirubin stains various shades of green, muscle and cell cytoplasm yellow, and collagen red.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

### STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse 	Running DI Water	--	1	--	
4	Immerse	Fouchet's Reagent	--	5	--	Once complete, <b>rinse in running DI water (1 minute)</b> and continue.
5	Immerse	Van Gieson's Stain	--	5	--	Without rinsing, continue to next step.
6	Rinse	95% Reagent Alcohol	--	--	30-60	
7	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
8	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
9	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 219.
2. With modifications by AMTS R&D Department, 1979-2019.


**AUTOSTAINER CONFIGURATION AND NOTES:**

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 4 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 15 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	500	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	16	480	30	200ml Bath Autostainer	2	320	160
40ml Hellendahl Jar	12	480	40	400ml Bath Autostainer	1	340	340

**CE MARKINGS AND DESIGNATIONS:**

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By	<b>EC REP</b>	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	<b>CE</b>		
	Corrosive		Health Hazard		Environmentally Damaging	
GHS05		GHS08		GHS09		

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
 1330 Thurman St.  
 Lodi, CA 95240 USA  
 Tel 800 860 4073  
 Fax 209 368 4136  
 www.americanmastertech.com

**StatLab**  
 2090 Commerce Drive  
 McKinney, TX 75069 USA  
 Tel 800 442 3573  
 Fax 972 436 1369  
 www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 219.
2. With modifications by AMTS R&D Department, 1979-2019.

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**COMPOSANTS INCLUS:** Fouchet's Reagent, Van Gieson's Stain

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