Human Brain

## Heidenhain's Aniline Blue-Azan Stain Kit Procedure

**100ml Kit Item #:** KTHEI

**Liter Kit Item#:** N/A

**Pint Kit Item #:** N/A

**Gallon Kit Item#:** N/A

**Control Slide(s)**
**Item#**

Bone Marrow CSB0325P

Skeletal Muscle CSS0725P

Uterus CSU0325P

**Included Components**

Azocarmine B Solution 5% Phosphotungstic Acid

0.1% Aniline Alcohol Solution Aniline Blue, Heidenhain's

1% Acetic Alcohol Solution


**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify collagen, muscle, reticulum, cytoplasm, chromatin, osteocytes, and neuroglia. Chromatin, neuroglia and osteocytes stain red, collagen and reticulum blue, cytoplasm pink to blue, and muscle red to yellow.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

**STAINING PROCEDURE:** Preheating Required. See step 4 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse 	Running DI Water	--	1	--	
4	Waterbath	Filtered Azocarmine B Solution	56°	15	--	Immerse into preheated solution. Once complete, discard solution and <b>rinse in running DI water (1 minute)</b> and continue.
5	Immerse	0.1% Aniline Alcohol Solution	--	--	10-15	Differentiate until cytoplasm and connective tissue are pale pink and nuclei are well defined (use microscope).
6	Rinse	1% Acetic Alcohol Solution	--	--	10-15	To stop differentiation.
7	Immerse	5% Phosphotungstic Acid	--	10-15	--	Once complete, <b>rinse in running DI water (1 minute)</b> and continue.
8	Immerse	Aniline Blue, Heidenhain's	--	10-15	--	Until the finest connective tissue fibers are well defined (use microscope). Once complete, <b>rinse in running DI water (1 minute)</b> and continue.
9	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB; Theory and practice of Histotechnology; 1980; 190.
2. With modifications by AMTS R&D Department, 1979-2019.

**AUTOSTAINER CONFIGURATION AND NOTES:**

This stain kit is not available in the pint or larger size.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	10	132	14	250ml Glass Stain Dish	n/a	n/a	n/a
30ml Glass Coplin Jar	6	119	20	200ml Bath Autostainer	n/a	n/a	n/a
40ml Hellendahl Jar	5	132	27	400ml Bath Autostainer	n/a	n/a	n/a

**CE MARKINGS AND DESIGNATIONS:**

	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
	Batch Code		Use By		Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use			
	Flammable		Health Hazard			

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

**StatLab**  
2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Sheehan DC Hrapchak BB; Theory and practice of Histotechnology; 1980; 190.
2. With modifications by AMTS R&D Department, 1979-2019.

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
**COMPOSANTS INCLUS:** Azocarmine B Solution, 0.1% Aniline Alcohol Solution, 1% Acetic Alcohol Solution, 5% Phosphotungstic Acid, Aniline Blue, Heidenhain's

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**COMPONENTES INCLUIDOS:** Azocarmine B Solution, 0.1% Aniline Alcohol Solution, 1% Acetic Alcohol Solution, 5% Phosphotungstic Acid, Aniline Blue Heidenhain's

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