Colon Adenocarcinoma

Illuminate™ Harris Stain Kit Procedure

Illuminate™ Harris

Liter Kit Item#: KTILLHARLT

Gallon Kit Item#: KTILLHARGAL

Illuminate™ Harris Phloxinated

Liter Kit Item#: KTILLHARPBLT

Gallon Kit Item#: KTILLHARPBGAL

Control Slide(s)	Item#
Heart	CSU0325P
Skeletal Muscle	CSS0725P
Uterus	CSU0325P

Included Components

Illuminate™ Harris Hematoxylin	Illuminate™ Bluing
Illuminate™ Differentiating	Illuminate™ Eosin or Eosin Phloxinated

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify nuclei and cytoplasm. Nuclei are regressively stained blue to purple and cytoplasm, collagen, and muscle pink to red.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE:

Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Illuminate™ Harris Hematoxylin	--	5	--	Adjust time to desired intensity. Once complete, rinse in running tap water (1 minute) and continue.
5	Immerse	Illuminate™ Differentiating	--	--	30	Once complete, rinse in running tap water (1 minute) and continue.
6	Immerse	Illuminate™ Bluing	--	1	--	Once complete, rinse in running tap water (1 minute) and continue.
7	Immerse	70% Reagent Alcohol	--	1	--	
8	Immerse	Illuminate™ Eosin or Eosin Phloxinated	--	1	--	Adjust time to desired intensity.
9	Immerse	95% Reagent Alcohol	--	1	--	
10	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 250 - 251.
2. With modifications by AMTS R&D Department, 1979-2018.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 7 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 18 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	500	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	16	480	30	200ml Bath Autostainer	2	320	160
40ml Hellendahl Jar	12	480	40	400ml Bath Autostainer	1	340	340

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Health Hazard	CE

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 250 - 251.
2. With modifications by AMTS R&D Department, 1979-2018.

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
COMPOSANTS INCLUS: Illuminate™ Harris Hematoxylin, Illuminate™ Differentiating, Illuminate™ Bluing, Illuminate™ Eosin ou Eosin Phloxinated

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COMPONENTES INCLUIDOS: Illuminate™ Harris Hematoxylin, Illuminate™ Differentiating, Illuminate™ Bluing, Illuminate™ Eosin o Eosin Phloxinated

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