Human Kidney

Jones Basement Membrane & Reticulum Stain Kit Procedure

100ml Kit Item #: KTJON

Liter Kit Item#: N/A

Pint Kit Item #: KTJONPT

Gallon Kit Item#: N/A

Control Slide(s)

Item#

Included Components

Kidney (Normal) CSK0125P

0.25M Sodium Borate 0.5% Periodic Acid

Reticulum Fibers CSR0125P

0.2M Boric Acid 0.2% Gold Chloride

3% Methenamine 3% Sodium Thiosulfate

5% Silver Nitrate Nuclear Fast Red Stain

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify basement membrane and reticulum fibers. Basement membranes and reticulum fibers stain black and nuclei and background pink to red.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 2-3µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

WORKING METHENAMINE-SILVER SOLUTION PREPARATION: Prepare at time of use. Solution expires in one hour.

#	Action	Amount	Chemical/Reagent	Details
1	Add	4ml	0.25M Sodium Borate Solution	Into a chemically clean container or new/unused plasticware.
2	Add	6ml	0.2M Boric Acid Solution	Mix thoroughly.
3	Add	42ml	3% Methenamine Solution	Mix thoroughly.
4	Add	3ml	5% Silver Nitrate	Mix thoroughly.
5	Filter	--	Medium Filter Paper	Filter into a chemically clean container or new/unused plasticware.


STAINING PROCEDURE: Preheating Required. See step 5 for more information. Microwave Procedure located on opposite page.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	0.5% Periodic Acid	--	10	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Waterbath	Working Methenamine - Silver Solution	65-70°	20-25	--	Immerse slide into a jar of room temperature Working Methenamine-Silver Solution and place jar into a preheated waterbath. Develop until basement membrane and reticulum fibers are stained intensely black. Use a quick DI rinse and microscope to check for correct differentiation. Do not allow tissue to cool or uneven staining will occur. Once complete, rinse in running DI water (1 minute) and continue at step 6.
6	Immerse	0.2% Gold Chloride	--	1	--	Section should be gray-black. Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	3% Sodium Thiosulfate	--	1-2	--	Once complete, rinse in running DI water (1 minute) and continue.
8	Immerse	Nuclear Fast Red Stain	--	7	--	Once complete, rinse in running DI water (1 minute) and continue.
9	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media				

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; Second Edition; 186 -187, 1980.
2. Luna LG: Histopathologic Methods and Color Atlas of Special Stain and Tissue Artifacts; American Histolabs, Inc.; 400 - 402, 1992.
3. With modifications by AMTS R&D Department, 1979-2019.

MICROWAVE STAINING PROCEDURE:

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohols, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse 	Running DI Water	--	1	--	
4	Microwave	0.5% Periodic Acid	70-80°	--	--	Immerse slide in 0.5% periodic acid, heat, and then incubate for 5 minutes. Once complete, rinse in running DI water (1 minute) and continue.
5	Microwave	Working Methenamine – Silver Solution	70-80°	--	--	Immerse slide into a loosely capped jar of Working Methenamine Silver Solution and heat (do not boil). Incubate for 1-5 minutes until basement membrane and reticulum fibers are stained intensely black. Use a quick DI rinse and microscope to check for correct differentiation. Do not allow tissue to cool or uneven staining will occur. Once complete, rinse in running DI water (1 minute) and continue.
6	Immerse	0.2% Gold Chloride	--	1	--	Section should be gray-black. Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	3% Sodium Thiosulfate	--	1-2	--	Once complete, rinse in running DI water (1 minute) and continue.
8	Immerse	Nuclear Fast Red Stain	--	7	--	Once complete, rinse in running DI water (1 minute) and continue.
9	Dehydrate	Absolute Alcohols, 3 changes	--	1	--	Fresh changes. 1 minute each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media				

AUTOSTAINER CONFIGURATION AND NOTES:

This kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Methenamine Silver Solution preparation procedure yields 55ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 9 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 20 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	30	120	4	250ml Glass Stain Dish	2	60	30
30ml Glass Coplin Jar	20	160	8	200ml Bath Autostainer	3	90	30
40ml Hellendahl Jar	15	240	16	400ml Bath Autostainer	1	50	50

CE MARKINGS AND DESIGNATIONS:

	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
	Batch Code		Use By		Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use			
	Irritant		Health Hazard		Environmentally Damaging	

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St. • Lodi, CA 95240 USA
Tel 800 860 4073 • Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive • McKinney, TX 75069 USA
Tel 800 442 3573 • Fax 972 436 1369
www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; Second Edition; 186 -187, 1980.
2. Luna LG: Histopathologic Methods and Color Atlas of Special Stain and Tissue Artifacts; American Histolabs, Inc.; 400 - 402, 1992.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 0.25M Sodium Borate, 0.2M Boric Acid, 3% Methenamine, 5% Silver Nitrate, 0.5% Periodic Acid, 0.2% Gold Chloride, 3% Sodium Thiosulfate, Nuclear Fast Red Stain

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