Human Brain

Luxol Fast Blue – Cresyl Echt Violet Stain Kit Procedure

100ml Kit Item #: KTLFB	Liter Kit Item#: N/A
Pint Kit Item #: KTLFBPT	Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components	
Cerebellum	CSC0825P	Luxol Fast Blue Stain	70% Reagent Alcohol
		0.05% Lithium Carbonate	Cresyl Echt Violet Stain


PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify myelin sheath and Nissl substance. Myelin fibers stain blue to turquoise, Nissl substance purple to violet, and nuclei purple to violet.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 8-10µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE: Preheating Required. See step 4 for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Waterbath	Luxol Fast Blue Stain	60°	60	--	Immerse into preheated solution. Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	0.05% Lithium Carbonate	--	--	0.5	3-5, quick dips, for half a second each to differentiate gray and white matter.
6	Immerse	70% Reagent Alcohol, 2 changes	--	--	0.5	Half second dips, each change until the gray and white matter is clearly distinguished; nuclei should be colorless and myelin should be turquoise against a pale blue-gray background. Use a quick DI water rinse and microscope to check for correct differentiation. Once complete, rinse in running DI water (1 minute) and continue.
7	Immerse	Cresyl Echt Violet Stain	--	10	--	Once complete, rinse in running DI water (5-10 secs) and continue.
8	Immerse	70% Reagent Alcohol	--	--	0.5	3-5, quick dips, for half a second each.
9	Dehydrate	Absolute Alcohol, 3 changes	--	--	3	3 seconds each change with agitation (excess time will decolorize Cresyl Echt Violet Stain).
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

1. Kluver H Barrera EA: A Method for the combined staining of cells and fibers in the nervous system; J Neuropathol Exp Neurol, 1953, 12: 400-403.
2. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 252 - 253, 263 - 264.
3. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 7 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 18 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	225	9	250ml Glass Stain Dish	2	169	85
30ml Glass Coplin Jar	16	216	14	200ml Bath Autostainer	2	144	72
40ml Hellendahl Jar	12	216	18	400ml Bath Autostainer	1	153	153

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Flammable		Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Kluver H Barrera EA: A Method for the combined staining of cells and fibers in the nervous system; J Neuropathol Exp Neurol, 1953, 12: 400-403.
2. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 252 - 253, 263 - 264.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Luxol Fast Blue Stain, 0.05% Lithium Carbonate, 70% Reagent Alcohol, Cresyl Echt Violet Stain.

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