Human Liver

## Copper Stain Kit Procedure, Microwave

100ml Kit Item #: KTMCO

Liter Kit Item#: KTMCOLT

Pint Kit Item #: KTMCOPT

Gallon Kit Item#: N/A

Control Slide(s)

Item#

Included Components

Copper

CSC0325P

Acetate Buffer

Rhodanine Stock

Modified Mayer's Hematoxylin

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify copper. Copper stains light brown to red, cytoplasmic granules and nuclei blue.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

**WORKING RHODANINE SOLUTION PREPARATION:** Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	46ml	Acetate Buffer	Into a coplin jar or container.
2	Add	4ml	Rhodanine Stock	Agitate Rhodanine Stock thoroughly before adding to Acetate Buffer.

**STAINING PROCEDURE:** Preheating required for waterbath procedure. See step 4a for more information.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	For routine waterbath procedure follow step 4a. For microwave procedure, follow step 4b.					

4a	Waterbath	Working Rhodanine Solution	37°	1080	--	Immerse into preheated solution for 18 hours. Once complete, continue at step 5.
----	-----------	----------------------------	-----	------	----	--

4b	Microwave	Working Rhodanine Solution	80-90°	--	--	Immerse slide into a loosely capped jar of Working Rhodanine Solution and heat (do not boil). Remove from microwave, cap jar tightly, and allow solution to cool to room temperature while agitating jar periodically. After cooling, examine using a microscope; repeat heating/cooling cycles as needed to reach desired intensity. Once complete, continue at step 5.
----	-----------	----------------------------	--------	----	----	--

5	Rinse	 Acetate Buffer, 2 changes	--	1	--	1 minute each change.
6	Immerse	Modified Mayer's Hematoxylin	--	--	5-10	With agitation. Once complete, <b>rinse in Acetate Buffer (2 changes at 1 minute each)</b> and continue.
7	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
8	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
9	Coverslip	Permanent Mounting Media	--	--	--	

- Lindquist RR: Studies on the Pathogenesis of Hepatolenticular II: Cytochemical methods for the location of copper; Arch Pathol; 87; 1969, 370.
- With modifications by AMTS R&D Department, 1979-2019.

### AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Working Rhodanine Solution preparation procedure yields approximately 50ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 5 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 16 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	500	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	16	480	30	200ml Bath Autostainer	2	320	160
40ml Hellendahl Jar	12	480	40	400ml Bath Autostainer	1	340	340

### CE MARKINGS AND DESIGNATIONS:

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By	<b>EC REP</b>	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	<b>CE</b>		
	Flammable		Health Hazard			
GHS02		GHS08				

### CONTACT INFORMATION:

**American MasterTech Scientific**  
1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

**StatLab**  
2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Lindquist RR: Studies on the Pathogenesis of Hepatolenticular II: Cytochemical methods for the location of copper; Arch Pathol; 87; 1969, 370.
2. With modifications by AMTS R&D Department, 1979-2019.

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**COMPOSANTS INCLUS:** Acetate Buffer, Rhodanine Stock, Modified Mayer's Hematoxylin.

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