Human Lung Tissue

## McDonalds Gram™ Stain Kit Procedure

100ml Kit Item #: KTMGS

Liter Kit Item#: KTMGSLT

Pint Kit Item #: KTMGSPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#
Gram +/- Gram -	CSG0225P

Included Components	
Gentian Violet	Carbol Fuchsin
Universal Iodine Solution™	Optimized Tartrazine Counterstain
Gram's Decolorizer	

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify gram positive and gram negative bacteria. Gram positive bacteria stain blue, gram negative bacteria red, other tissue yellow, and nuclei red.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

### STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	Gentian Violet	--	1	--	Once complete, <b>rinse in running DI water (1 minute)</b> to remove excess dye and continue.
5	Immerse	Universal Iodine Solution™	--	1	--	Once complete, <b>rinse in running DI water (1 minute)</b> to remove excess iodine and continue.
6	Decolorize	Gram's Decolorizer	--	--	5-15	Until color does not bleed from tissue. Once complete, <b>rinse in running DI water (30 seconds)</b> and continue.
7	Immerse	Carbol Fuchsin	--	--	90	Once complete, <b>rinse in running DI water (15 seconds)</b> and continue.
8	Immerse	Optimized Tartrazine Counterstain	--	1	--	
9	Dehydrate	Absolute Alcohol, 3 changes	--	--	3	3 seconds each change with agitation.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 235
2. With modifications by AMTS R&D Department, 1979-2019.

**AUTOSTAINER CONFIGURATION AND NOTES:**

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 6 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 17 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	225	9	250ml Glass Stain Dish	2	169	85
30ml Glass Coplin Jar	16	216	14	200ml Bath Autostainer	2	144	72
40ml Hellendahl Jar	12	216	18	400ml Bath Autostainer	1	153	153

**CE MARKINGS AND DESIGNATIONS:**

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By	<b>EC REP</b>	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	<b>CE</b>		
	Flammable		Irritant		Health Hazard	

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

**StatLab**  
2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 235
2. With modifications by AMTS R&D Department, 1979-2019.

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**COMPOSANTS INCLUS:** Gentian Violet, Universal Iodine Solution™, Gram's Decolorizer, Carbol Fuchsin, Optimized Tartrazine Counterstain

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**COMPONENTES INCLUIDOS:** Gentian Violet, Universal Iodine Solution™, Gram's Decolorizer, Carbol Fuchsin, Optimized Tartrazine Counterstain

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