Human Bone Marrow

Optimized Giemsa™ Stain Kit Procedure

100ml Kit Item #: N/A	Liter Kit Item#: N/A
Pint Kit Item #: KTOGIPT	Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components
Bone Marrow	CSB0325P	Optimized Giemsa Buffer™
Helicobacter Pylori	CSH0125P	Giemsa Stain (Stock)
Mast Cell	CSM0125P	
Toxoplasma Gondii	CST0325P	

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify bone marrow cells and helicobacter. Nuclei is stained blue to violet, rickettsias intense reddish-purple, cytoplasm light blue, collagen and muscle pale pink, erythrocytes gray, yellow or pink, helicobacter blue, and toxoplasma gondii blue.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

WORKING GIEMSA STAIN PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	50ml	Optimized Giemsa Buffer™	Into a chemically cleaned container or new/unused plasticware.
2	Add	2.5ml	Giemsa Stain (Stock)	Mix thoroughly.

STAINING PROCEDURE:

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Optimized Giemsa Buffer™	--	3	--	
5	Immerse	Working Giemsa Stain	--	7-11	--	Agitate every few minutes to ensure proper staining.
6	Rinse	Optimized Giemsa Buffer™	--	--	10	Until any remaining stain runs off.
7	Immerse	Optimized Giemsa Buffer™	--	3	--	
8	Air Dry	On Benchtop, or in Empty Bath	--	5-10	--	Until section is completely dry.
9	Clear	Xylene or Substitute	--	--	1	5-10 dips for 1 second each.
10	Coverslip	Permanent Mounting Media	--	--	--	

1. A.F.I.P. Manual: pg 212. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 155-157
2. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Working Giemsa Stain preparation procedure yields approximately 53ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 5 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 12 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	31	620	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	20	600	30	200ml Bath Autostainer	3	480	160
40ml Hellendahl Jar	15	600	40	400ml Bath Autostainer	1	340	340

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Flammable		Toxic		Health Hazard	
GHS02		GHS06		GHS08		

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. A.F.I.P. Manual: pg 212. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 155-157
2. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Optimized Giemsa Buffer™, Giemsa Stain (Stock)

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