Human Cervical Smear

PAP Stain Kit Procedure, MasterTech Modified

100ml Kit Item #: N/A	Liter Kit Item#: KTPAPLT
Pint Kit Item #: KTPAPPT	Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components
N/A		Modified Mayer's Hematoxylin OG-6 Cytology Stain EA-50 Cytology Stain


PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared cytology smears (in vitro) to identify atypical, inflammatory, and malignant cells in cytology smears. Chromatin stains blue, keratin orange, superficial squamous cells various shades of pink, nucleoli, cilia, red blood cells various shades of pink, and all metabolic cell cytoplasm various shades of blue-green.

SPECIMEN CRITERIA: A cytology smear fixed in 95% Reagent Alcohol or cytology fixative.

STORAGE AND USAGE NOTES: If Harris or Gill's Hematoxylin is used, add differentiating and bluing steps. Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Hydrate	70% Reagent Alcohol, 2 changes	--	5	--	5 minutes each change.
2	Hydrate	50% Reagent Alcohol, 1 change	--	1	--	1 minute each change.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	Modified Mayer's Hematoxylin	--	3-5	--	
5	Rinse	Running Tap Water	--	2	--	
6	Immerse 	95% Reagent Alcohol, 2 changes	--	--	20	20 seconds each change with gentle agitation.
7	Immerse	OG-6 Cytology Stain	--	2	--	With gentle agitation. Once complete, immerse in 95% Reagent Alcohol (2 changes at 15 seconds each) and continue.
8	Immerse	EA-50 Cytology Stain	--	3	--	With gentle agitation.
9	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
10	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
11	Coverslip	Permanent Mounting Media	--	--	--	

1. Carson FL, Hladik C Histotechnology: A Self-Instructional Text (3rd ed.)
2. American Society for Clinical Pathology Press [2009]. 362-363.
3. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 7 baths is required to perform this procedure excluding hydration, dehydration, and clearing, or 16 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	500	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	16	480	30	200ml Bath Autostainer	2	320	160
40ml Hellendahl Jar	12	480	40	400ml Bath Autostainer	1	340	340

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Flammable		Health Hazard			

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Carson FL, Hladik C Histotechnology: A Self-Instructional Text (3rd ed.)
2. American Society for Clinical Pathology Press [2009]. 362-363.
3. With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: Modified Mayer's Hematoxylin, OG-6 Cytology Stain, EA-50 Cytology Stain

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