Human Colon

## Periodic Acid-Schiff's (P.A.S.) Stain Kit Procedure

**100ml Kit Item #:** KTPAS

**Liter Kit Item#:** KTPASLT

**Pint Kit Item #:** KTPASPT

**Gallon Kit Item#:** KTPASGAL

| Control Slide(s) | Item#    |
|------------------|----------|
| Aspergillus      | CSA0625P |
| Candida Albicans | CSC0525P |
| Glycogen         | CSG0125P |

### Included Components

|                             |                              |
|-----------------------------|------------------------------|
| 0.5% Periodic Acid          | Modified Mayer's Hematoxylin |
| Optimized Schiff's Solution | Light Green Counterstain     |

### Components Not Included:

Diastase (Only required for Glycogen digestion) Item# AHDIA50

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify fungi, glycogen and mucin. Basement membrane, fungi, glycogen and mucin stain pink to red, other tissue green, and nuclei blue if hematoxylin is used.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

**STORAGE AND USAGE NOTES:** (Diastase, AMTS Item# AHDIA50) is available from American MasterTech. Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

**DIASTASE WORKING SOLUTION:** Prepare only if digestion is desired! Solution expires after one use.

| # | Action | Amount | Chemical/Reagent                 | Details  |
|---|--------|--------|----------------------------------|--|
| 1 | Add    | 1gm    | Diastase ( <i>Not Included</i> ) | Into a chemically cleaned container or new/unused plasticware. |
| 2 | Add    | 100ml  | Distilled/DI Water               | Mix thoroughly.  |

### STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

| # | Action   | With                            | Heat °C | Time |      | Details  |
|---|--|---------------------------------|---------|------|------|--|
|   |  |                                 |         | Mins | Secs |  |
| 1 | Deparaffinize  | Xylene or Substitute, 2 changes | --      | 5    | --   | 5 minutes each change or as required if using a xylene substitute. |
| 2 | Rinse  | Absolute Alcohol, 3 changes     | --      | 1    | --   | 1 minute each change or as required if using graded alcohols.      |
| 3 | Rinse  | Running DI Water                | --      | 1    | --   |  |
| 4 | If digestion is desired, follow step 4a below. If not, continue to step 5. |                                 |         |      |      |  |

#### 4a. Digestion Step

|         |                           |    |    |    |  |
|---------|---------------------------|----|----|----|--|
| Immerse | Diastase Working Solution | -- | 20 | -- | Once complete, <b>rinse in running DI water (1 minute)</b> and continue. |
|---------|---------------------------|----|----|----|--|

|    |           |                                 |    |    |       |  |
|----|-----------|---------------------------------|----|----|-------|--|
| 5  | Immerse   | 0.5% Periodic Acid              | -- | 7  | --    | Once complete, <b>rinse in running DI water (10 seconds)</b> and continue. |
| 6  | Immerse   | Optimized Schiff's Solution     | -- | 15 | --    |  |
| 7  | Rinse     | Running Tap Water               | -- | 5  | --    |  |
| 8  | Immerse   | Modified Mayer's Hematoxylin    | -- | 2  | --    | Once complete, <b>rinse in running tap water (1 minute)</b> and continue.  |
| 9  | Immerse   | Light Green Counterstain        | -- | -- | 30-60 | Until desired background intensity is reached.                             |
| 10 | Dehydrate | Absolute Alcohol, 3 changes     | -- | 1  | --    | 1 minute each change.  |
| 11 | Clear     | Xylene or Substitute, 3 changes | -- | 1  | --    | 1 minute each change or as required if using a xylene substitute.          |
| 12 | Coverslip | Permanent Mounting Media        | -- | -- | --    |  |

1. Mac Manus: Staining Technology; pg 99.
2. With modifications by AMTS R&D Department, 1979-2019.

**AUTOSTAINER CONFIGURATION AND NOTES:**

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 6 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, clearing, digestion, and nuclear staining or 18 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

| Bath Type              | Uses | Slides | S.C. | Bath Type              | Uses | Slides | S.C. |
|------------------------|------|--------|------|------------------------|------|--------|------|
| 20ml Plastic Slide Jar | 25   | 500    | 20   | 250ml Glass Stain Dish | 2    | 375    | 188  |
| 30ml Glass Coplin Jar  | 16   | 480    | 30   | 200ml Bath Autostainer | 2    | 320    | 160  |
| 40ml Hellendahl Jar    | 12   | 480    | 40   | 400ml Bath Autostainer | 1    | 340    | 340  |

**CE MARKINGS AND DESIGNATIONS:**

|  |                                    |  |                                   |   |                |  |
|--|------------------------------------|--|-----------------------------------|---|----------------|--|
| <b>REF</b>   | Catalogue Number                   |   | Temperature Limitation            |  | Manufacturer   | American MasterTech Scientific<br>1330 Thurman St.<br>Lodi, CA 95240 USA<br>Tel 800 860 4073<br>Fax 209 368 4136 |
| <b>LOT</b>   | Batch Code                         |   | Use By                            | <b>EC REP</b>   | Representative | Emergo Europe<br>Prinsessegracht 20<br>2514 AP The Hague<br>The Netherlands                                      |
| <b>IVD</b>   | In Vitro Diagnostic Medical Device |   | Consult Instructions Prior to Use | <b>CE</b>   |                |  |
|  | Irritant                           |  | Health Hazard                     |   |                |  |

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

**StatLab**  
2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Mac Manus: Staining Technology; pg 99.
2. With modifications by AMTS R&D Department, 1979-2019.

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