Human Lung

Pneumocystis CN 1.5 Stain Kit Procedure

100ml Kit Item#: KTPNE

Liter Kit Item#: N/A

Pint Kit Item#: KTPNEPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#
Pneumocystis	CSP0125P

Included Components	
Cresyl Echt Violet	Slide Jar(s)
Naphthol Yellow S	

Components Not Included:
Glacial Acetic Acid (Item#: SPA0243)
Sulfuric Acid (Item#: SPS2743)

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify *Pneumocystis jirovecii* organisms. *Pneumocystis jirovecii* stain violet to purple, connective tissue blue to green, erythrocytes yellow, and mucin and cartilage rose to purple.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 4-5µm tissue sections.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

WORKING Sulfation REAGENT PREPARATION: Prepare solution at time of use. Wear protective clothing, gloves, and goggles while making this solution. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	15ml	Glacial Acetic Acid (<i>Not Included</i>)	Into slide jar provided in this kit.
2	Add	5ml	Sulfuric Acid (<i>Not Included</i>)	Very slowly add to slide jar. Place cap on tightly and invert slide jar several times to mix the acids. Wait 5-10 minutes for the heat generated during mixing to dissipate before placing slide into jar. Use the slide jar for Sulfation Reagent only.

STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running Tap Water	--	1	--	
4	Immerse	Working Sulfation Reagent	--	10	--	Agitate periodically to keep Acetic Acid from settling or Pneumocystis will stain poorly. Do not stain more than two slides at a time. Once complete, rinse in running tap water (1 minute) and continue.
5	Immerse	Cresyl Echt Violet	--	10-15	--	Agitate periodically. Pneumocystis should be stained violet (use microscope). Once complete, rinse in running tap water (5 seconds) and continue.
6	Immerse	Naphthol Yellow S	--	--	4	Excessive counterstain time will decolorize the Pneumocystis.
7	Dehydrate	Absolute Alcohol	--	--	1	2-3 dips for 1 second each.
8	Air Dry	On Benchtop, or in Empty Bath	--	1-2	--	Until section is completely dry.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	Fresh changes. 1 minute each change.
10	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 250 - 251.
2. With modifications by AMTS R&D Department, 1979-2019.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit in the pint size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Sulfation Reagent preparation procedure yields approximately 20ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 5 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 10 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	500	20	250ml Glass Stain Dish	2	375	188
30ml Glass Coplin Jar	16	480	30	200ml Bath Autostainer	2	320	160
40ml Hellendahl Jar	12	480	40	400ml Bath Autostainer	1	340	340

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Health Hazard	CE

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 250 - 251.
2. With modifications by AMTS R&D Department, 1979-2019.

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