Animal Testes

Steiner Stain Kit Procedure

100ml Kit Item #: KTSTE

Liter Kit Item#: KTSTELT

Pint Kit Item #: KTSTEPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components	
Aspergillus	CSA0625P	Hydroquinone Capsules	Oxidizer Solution
Candida Albicans	CSC0525P	2.5% Gum Mastic	10% Zinc Formalin
Cat Scratch Fever	CSC0225P	1% Silver Nitrate	0.1% Silver Nitrate
Helicobacter pylori	CSH0125P		
Legionnaires	CSL0125P		
Spirochete	CSS0125P		

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify fungi, Helicobacter pylori, Cat scratch fever, Legionella pneumophila, and spirochetes. Cat scratch fever, Legionella pneumophila, spirochetes, Helicobacter pylori and fungi stain dark brown, gray or black, and the background stains yellow, tan or brown.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 3-4µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

REDUCING SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	25ml	Distilled Water	Into a chemically cleaned container or new/unused plasticware.
2	Add	1 capsule	Hydroquinone Capsule	Empty contents of capsule and mix until dissolved.
3	Add	15ml	2.5% Gum Mastic	Mix thoroughly.
4	Filter	--	Medium Filter Paper	Filter into a chemically cleaned container or new/unused plasticware.
5	Add	6 Drops	1% Silver Nitrate	Mix thoroughly.

STAINING PROCEDURE: Preheating Required. See steps 6 and 10 for more information. Microwave procedure located on opposite page.

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	Oxidizer Solution	--	20	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	10% Zinc Formalin	--	5	--	Once complete, rinse in running DI water (1 minute) and continue.
6	Waterbath	0.1% Silver Nitrate	70°	5	--	Immerse into preheated solution. Once complete, discard solution and rinse in running DI water (1 minute) and continue.
7	Dehydrate	Absolute Alcohol, 2 changes	--	1	--	1 minute each change with agitation.
8	Immerse	2.5% Gum Mastic	--	3	--	Without rinsing, continue to next step.
9	Air Dry	On Benchtop, or in Empty Bath	--	1-2	--	Until gum mastic is completely dry.
10	Waterbath	Reducing Solution	70°	10-15	--	Immerse into preheated solution until section is tan to brown. Once complete, discard solution.
11	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
12	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
13	Coverslip	Permanent Mounting Media	--	--	--	

- Elias JM, Green C: Modified Steiner method for the demonstration of Spirochetes in tissue, Am J Clin Pathol; 1979, 71: 109 - 111.
- Margeson LS, Chapman CM: Use of Zinc Formalin as a Sensitizer in Silver Stain For Spirochetes; The Journal of Histotechnology; Vol. 19, No. 2; June 1996, 135-138.
- With modifications by AMTS R&D Department, 1979-2018.

MICROWAVE STAINING PROCEDURE:

 Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohols, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	Oxidizer Solution	--	20	--	Once complete, rinse in running DI water (1 minute) and continue.
5	Immerse	10% Zinc Formalin	--	5	--	Once complete, rinse in running DI water (1 minute) and continue.
6	Microwave	0.1% Silver Nitrate	80-90°	--	--	Heat 0.1% Silver Nitrate in a loosely capped slide jar (without slide) Do not boil. Remove from microwave, tighten cap, and agitate to equalize temperature. Immerse slide, agitate thoroughly, and incubate for 2 minutes. Once complete, rinse in running DI water (1 minute) and continue.
7	Dehydrate	Absolute Alcohol, 2 changes	--	1	--	1 minute each change.
8	Immerse	2.5% Gum Mastic	--	3	--	During this step, prepare the Reducing Solution.
9	Air Dry	On Benchtop, or in Empty Bath	--	1-2	--	Until gum mastic is completely dry.
10	Microwave	Reducing Solution	80-90°	--	--	Heat Reducing Solution in a loosely capped slide (without slide). Do not boil. Remove from microwave, tighten cap, and agitate to equalize temperature. Immerse slide, agitate thoroughly, and incubate until section is brown.
11	Dehydrate	Absolute Alcohols, 3 changes	--	1	--	Fresh changes. 1 minute each change.
12	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
13	Coverslip	Permanent Mounting Media	--	--	--	Coverslip immediately.

AUTOSTAINER CONFIGURATION AND NOTES:

This kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. Standard Reducing Solution preparation procedure yields 40ml of solution and must be scaled up to accommodate desired autostainer bath size. A minimum of 8 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 20 baths to run the complete procedure.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	100	4	250ml Glass Stain Dish	2	60	30
30ml Glass Coplin Jar	16	128	8	200ml Bath Autostainer	2	60	30
40ml Hellendahl Jar	12	192	16	400ml Bath Autostainer	1	50	50

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	CE		
	Flammable		Corrosive		Irritant	
	Health Hazard		Environmentally Damaging			

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St. • Lodi, CA 95240 USA
Tel 800 860 4073 • Fax 209 368 4136 • www.americanmastertech.com

StatLab
2090 Commerce Drive • McKinney, TX 75069 USA
Tel 800 442 3573 • Fax 972 436 1369 • www.statlab.com

- Elias JM, Green C: Modified Steiner method for the demonstration of Spirochetes in tissue, Am J Clin Pathol; 1979, 71: 109 - 111.
- Margeson LS, Chapman CM: Use of Zinc Formalin as a Sensitizer in Silver Stain For Spirochetes; The Journal of Histotechnology; Vol. 19, No. 2; June 1996, 135-138.
- With modifications by AMTS R&D Department, 1979-2018.

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COMPOSANTS INCLUS: Hydroquinone Capsules, 2.5% Gum Mastic, 1% Silver Nitrate, Oxidizer Solution, 10% Zinc Formalin, 0.1% Silver Nitrate.

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1. Elias JM, Green C: Modified Steiner method for the demonstration of Spirochetes in tissue, Am J Clin Pathol; 1979, 71: 109 - 111.
2. Margeson LS, Chapman CM: Use of Zinc Formalin as a Sensitizer in Silver Stain For Spirochetes; The Journal of Histotechnology; Vol. 19, No. 2; June 1996, 135-138.
3. With modifications by AMTS R&D Department, 1979-2018.

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COMPONENTES INCLUIDOS: Hydroquinone Capsules, 2.5% Gum Mastic, 1% Silver Nitrate, Oxidizer Solution, 10% Zinc Formalin, 0.1% Silver Nitrate.

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