Human Kidney

## Von Kossa Stain Kit Procedure

**100ml Kit Item#:** KTVKO      **Liter Kit Item#:** KTVKOLT  
**Pint Kit Item#:** KTVKOPT      **Gallon Kit Item#:** N/A

Control Slide(s)	Item#	Included Components
Calcium	CSC0125P	5% Silver Nitrate 5% Sodium Thiosulfate Nuclear Fast Red Stain

**PRINCIPLE AND RESULTS:** This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify calcium deposits. Calcium salts stain black, cytoplasm pink, and nuclei red.

**SPECIMEN CRITERIA:** Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

**STORAGE AND USAGE NOTES:** Store/Use each component according to the temperature and expiration on the label.

**PRECAUTIONS:** For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

### STAINING PROCEDURE:

Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	Immerse	5% Silver Nitrate	--	10-60	--	Expose to an ultraviolet light (10-20 minutes) or 100 watt incandescent desk lamp (30-60 minutes). Once complete, <b>rinse in running DI water (1 minute)</b> and continue.
5	Immerse	5% Sodium Thiosulfate	--	2-3	--	Once complete, <b>rinse in running DI water (1 minute)</b> and continue.
6	Immerse	Nuclear Fast Red Stain	--	5	--	Once complete, <b>rinse in running DI water (1 minute)</b> and continue.
7	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	Fresh changes. 1 minute each change.
8	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
9	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 227.
2. With modifications by AMTS R&D Department, 1979-2019.

**AUTOSTAINER CONFIGURATION AND NOTES:**

This stain kit in the pint or larger size may be easily adapted for use on most open-platform autostainers using the staining procedure grid on the reverse side of this page. A minimum of 4 baths is required to perform this procedure excluding deparaffinization, hydration, dehydration, and clearing, or 15 baths to run the complete procedure.

**TEST YIELD:** \*Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	25	100	4	250ml Glass Stain Dish	2	60	30
30ml Glass Coplin Jar	16	128	8	200ml Bath Autostainer	2	60	30
40ml Hellendahl Jar	12	192	16	400ml Bath Autostainer	1	50	50

**CE MARKINGS AND DESIGNATIONS:**

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By	<b>EC REP</b>	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	<b>CE</b>		
	Irritant		Health Hazard		Environmentally Damaging	

**CONTACT INFORMATION:**

**American MasterTech Scientific**  
1330 Thurman St.  
Lodi, CA 95240 USA  
Tel 800 860 4073  
Fax 209 368 4136  
www.americanmastertech.com

**StatLab**  
2090 Commerce Drive  
McKinney, TX 75069 USA  
Tel 800 442 3573  
Fax 972 436 1369  
www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 227.
2. With modifications by AMTS R&D Department, 1979-2019.

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**COMPOSANTS INCLUS:** 5% Silver Nitrate, 5% Sodium Thiosulfate, Nuclear Fast Red Stain

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