Human Gastric Tissue

Modified Stat Warthin-Starry Stain Kit Procedure

100ml Kit Item #: KTWST

Liter Kit Item#: N/A

Pint Kit Item #: KTWSTPT

Gallon Kit Item#: N/A

Control Slide(s)	Item#	Included Components	
Cat Scratch Fever	CSC0225P	4% Gelatin Acidulated	Spirochete Sensitizer
Helicobacter pylori	CSH0125P	2% Silver Nitrate Acidulated	0.5% Silver Nitrate Acidulated
Legionnaires	CSL0125P	0.1% Hydroquinone Acidulated	
Spirochete	CSS0125P		

PRINCIPLE AND RESULTS: This kit is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify Helicobacter pylori, spirochetes, Legionella pneumophila, and Cat Scratch Fever bacteria. Helicobacter pylori, Legionella pneumophila, spirochetes, and Cat Scratch Fever bacteria stain dark brown, gray or black and background yellow to brown.

SPECIMEN CRITERIA AND: Appropriately fixed, paraffin-embedded, 4-5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

REDUCING SOLUTION PREPARATION: Prepare solution at time of use. Solution expires after one use.

#	Action	Amount	Chemical/Reagent	Details
1	Add	12.5ml	4% Gelatin Acidulated	Warm 4% Gelatin Acidulated to liquefy. Place into an unused container.
2	Add	20-30 Drops	2% Silver Nitrate Acidulated	Mix thoroughly.
3	Add	7.5ml	0.1% Hydroquinone Acidulated	Mix thoroughly or uneven staining will occur.

STAINING PROCEDURE: Preheating Required. See steps 5 and 6 for more information. Microwave procedure located on opposite page.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	For spirochete staining, follow step 4a and 4b below. For all other bacteria, continue to step 5.					

4a	Immerse	Spirochete Sensitizer	--	5	--	
4b	Rinse	Running DI Water	--	1	--	Once complete, continue at step 5.

5	Waterbath	0.5% Silver Nitrate	65-70°	3-5	--	Immerse into preheated solution. Agitate periodically. Without rinsing, continue to next step.
6	Waterbath	Reducing Solution	65-70°	1-5	--	Immerse slide into a jar of room temperature Reducing Solution and place jar into a preheated waterbath. Develop until section is brown. Agitate frequently (or uneven staining will occur).
7	Rinse	Hot Running Tap Water	--	2	--	Rinse carefully in hot tap water.
8	Dehydrate	Absolute Alcohol, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	Coverslip immediately.

1. Kerr DA: Improved Warthin-Starry method for tissue sections; Am J Clin Pathol; 1938, 8: 63 - 67.
2. With modifications by AMTS R&D Department, 1979-2019.

MICROWAVE STAINING PROCEDURE:

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohols, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	Running DI Water	--	1	--	
4	For spirochete staining, follow step 4a and 4b below. For all other bacteria, continue to step 5.					
4a	Immerse	Spirochete Sensitizer	--	5	--	
4b	Rinse	Running DI Water	--	1	--	Once complete, continue at step 5.
5	Microwave	0.5% Silver Nitrate	80-90°	--	--	Heat 0.5% Silver Nitrate in a loosely capped slide jar (without slides). Do not boil. Remove from microwave, tighten cap, and agitate to equalize temperature. Immerse slide, agitate thoroughly, and incubate for 5 minutes.
6	Microwave	Reducing Solution	80-90°	--	--	Heat Reducing Solution in a loosely capped slide jar (without slides). Do not boil. Transfer slide from 0.5% Silver Nitrate into the Reducing Solution and agitate thoroughly (or uneven staining will occur). Incubate for 2 to 5 minutes; remove slide when section is brown.
7	Rinse	Hot Running Tap Water	--	2	--	Rinse carefully in hot tap water.
8	Dehydrate	Absolute Alcohols, 3 changes	--	1	--	1 minute each change.
9	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
10	Coverslip	Permanent Mounting Media	--	--	--	Coverslip immediately.

AUTOSTAINER CONFIGURATION AND NOTES:

This stain kit is not adapted for use on open-platform autostainers.

TEST YIELD: *Assumes pint kit and maximum slides per run. Actual Results may vary. S.C. denotes number of slides between "Solution Change".

Bath Type	Uses	Slides	S.C.	Bath Type	Uses	Slides	S.C.
20ml Plastic Slide Jar	32	128	4	250ml Glass Stain Dish	2	60	30
30ml Glass Coplin Jar	21	168	8	200ml Bath Autostainer	n/a	n/a	n/a
40ml Hellendahl Jar	16	256	16	400ml Bath Autostainer	n/a	n/a	n/a

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By	EC REP	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
IVD	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use		Irritant	

CONTACT INFORMATION:

American MasterTech Scientific
 1330 Thurman St.
 Lodi, CA 95240 USA
 Tel 800 860 4073
 Fax 209 368 4136
 www.americanmastertech.com

StatLab
 2090 Commerce Drive
 McKinney, TX 75069 USA
 Tel 800 442 3573
 Fax 972 436 1369
 www.statlab.com

- Kerr DA: Improved Warthin-Starry method for tissue sections; Am J Clin Pathol; 1938, 8: 63 - 67.
- With modifications by AMTS R&D Department, 1979-2019.

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COMPOSANTS INCLUS: 4% Gelatin Acidulated, 2% Silver Nitrate Acidulated, 0.1% Hydroquinone Acidulated, Spirochete Sensitizer, 0.5% Silver Nitrate Acidulated

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