American MasterTech  
scientific laboratory supplies

# CoverSafe™

## Low-Hazard Mounting Medium

(Rev. 3/8/18)

### **CoverSafe™ The product and its purpose:**

American MasterTech Scientific, Inc. has been asked many times by Histologists, Cytologists, and Pathologists to develop a low-hazard mounting medium. We have developed a product without Xylene, Toluene, Dioxane, or Benzene as its base solvent. CoverSafe™ is the result of extensive research, experimentation, and advice from some of the most advanced polymer research institutions in the world. CoverSafe™ is unlike any other permanent mounting medium available today; it is low-hazard, fast-drying, and optically clear. CoverSafe™ works well with all clearing agents and because it is non-volatile you will not have problems with bubbles, solvent evaporation spaces, or mounting medium flow. CoverSafe™ dries quickly and retains excellent refractivity indefinitely. CoverSafe™ can be used in automatic coverslipping machines. CoverSafe™ does not contain chemicals listed in OSHA or Federal EPA hazardous materials regulations. CoverSafe™ is intended for In Vitro Diagnostic use. Active ingredients include; D-Limonene and Acrylic Copolymer.

### **CoverSafe™, Instructions for use:**

Use CoverSafe™ as you would any permanent mounting medium; for best results, apply an extra drop or two of CoverSafe™ and blot slide thoroughly. Drying time will vary depending on room temperature and humidity. To accelerate drying, place slides in a laboratory oven set at 60 - 75°C for 15 to 30 minutes. To remove a coverslip, soak slide in Xylene until coverslip detaches. If desired, CoverSafe™ viscosity may be reduced by adding a small amount of D-Limonene or Xylene. NOTE: Picric acid stains may bleed or fade.

### **CoverSafe™, Storage and Stability:**

Store CoverSafe™ at 15 - 30°C; it is stable under normal laboratory conditions. Keep tightly sealed when not in use. See product label for shelf life.

### **CoverSafe™, Precautions:**

See Safety Data Sheet (SDS) for warnings and precautions.

### **CoverSafe™, Disposal:**

Dispose of CoverSafe™ in accordance with all local, state, federal, environmental and/or EU regulations.

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**CE MARKINGS AND DESIGNATIONS:**

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By		Representative	Emergo Europe Prinsesegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	 GHS02	Flammable	
 GHS07	Irritant	 GHS09	Environmentally damaging			

**Copyright ©2015-2018 American MasterTech Scientific. All rights reserved.**

*No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.*