



MasterTech LYSING SOLUTION PROCEDURE

Item# CYC01

(Revised 03/8/18)

PRINCIPLE: Rapid lysing of red blood cells contained in Cytology and Histology specimens.

PROCEDURE:

1. Concentrate specimen by centrifuging as usual.
2. Carefully pour off supernatant as usual.
3. Add 5 ml of MasterTech Lysing Solution for every 1 ml of specimen.
4. Vortex solution to suspend cells.
5. Allow suspension to stand for 1 hour.
6. Centrifuge suspension to concentrate cells.
7. Pour off supernatant liquid.
8. Wrap cell block as usual and place in tissue cassette.
9. Place cassette 10% Formalin (or fixative of choice) and process as usual.

Copyright ©2007-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.