

# PolyCyte™ Mono-Layer Cytology Fixative Procedure

Item#: CYP CYCS, CYP CYLT

(Revised 03/08/18)



American MasterTech  
scientific laboratory supplies

**PRINCIPLE:** PolyCyte™ allows mono-layer cell preparation for cytology staining.

**SPECIMEN:** GYN, FNA, cell buttons, urine cytology, or other cytology specimens.

**QUALITY CONTROL:** American MasterTech Recommended Control Slide: Not Applicable.


**COLLECTION METHODS:** For GYN specimens, place collected cells into PolyCyte™ pre-filled bottle by agitating cytology-brush or cytology-broom in the container. For break-off brushes, agitate brush then break off and leave in PolyCyte™ container. For all other cellular specimens, follow standard collection protocol then deposit cellular material into the PolyCyte™ container.

**Important Note:** Insure that the cap on the PolyCyte™ container is properly affixed and tightened before sending the container to the laboratory.

**PROCEDURE:** (PolyCyte™ Spray Procedure located on opposite page.)

- 1 Vortex or spin pre-filled PolyCyte™ container for 5 to 10 seconds. For specimens containing scanty amount of cellular material see **procedure note** below.
- 2 Place SuperFrost Plus™ slides on a flat surface. Using a disposable pipet, dispense 3 to 4 drops of the PolyCyte™ suspended specimen at the same spot on the slides. Keeping slides flat, allow PolyCyte™ suspended specimen to dry at room temperature for 15 to 30 minutes. (Longer time at room temperature will improve mono-layer.)
- 3 Keeping slide flat, place slides in a 80°C pre-heated laboratory oven for 10 minutes. Be sure that slides are completely dry before removing them from the oven!
- 4 Remove slides from oven and allow them to cool until they reach room temperature.
- 5 Hydrate slides through 2 changes of 50% Alcohol for 3 to 5 minutes each with gentle agitation.
- 6 Place slides in stain rack.
- 7 Stain slides using Papanicolaou stain procedure.

**PROCEDURE NOTE:** For specimens containing a scanty amount of cellular material, vortex or spin pre-filled PolyCyte™ container for 5-10 seconds. Decant into a vacutainer using evacuation straw and centrifuge specimen for 1 to 5 minutes at 3000 RPM. Next, pour off supernatant, being careful to not pour off cells, add 2ml of fresh PolyCyte™ to the container, and vortex vacutainer again to re-mix the cells. Continue to step two.



Cytology Smear

Copyright © 2012-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.

# PolyCyte™ Mono-Layer Cytology Spray Procedure

Item#: CYPCY40Z

(Revised 03/08/18)

**PRINCIPLE:** PolyCyte™ allows mono-layer cell preparation for cytology staining.

**SPECIMEN:** GYN, FNA, cell buttons, urine cytology, or other cytology specimens.


**QUALITY CONTROL:** American MasterTech Recommended Control Slide: Not Applicable.

## PROCEDURE:

- 1 Collect sample with a cytology brush or broom.
- 2 Smear slide with sample.
- 3 Spray PolyCyte™ onto sample from a distance of 5cm.
- 4 Allow slide to dry at room temperature for 15 minutes.
- 5 Hydrate slide through 2 changes of 50% Alcohol for 3 minutes with gentle agitation.
- 6 Place slide on stain rack.
- 7 Stain slide with stain procedure.



American MasterTech  
scientific laboratory supplies



Cytology Smear

Copyright © 2012-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.