Living up to Life

Surgipath®

Paraplast®

Paraplast® Plus™

Paraplast® X-tra™

Blue Ribbon™

EM-400® Embedding Medium

Paraplast® Bulk

Paraplast® High Melt

Formula "R"™

Tissue Infiltration Medium

Tissue Embedding Media

English.....	2
Български (Bulgarian).....	3
中文 (简体) (Chinese).....	4
Čeština (Czech).....	5
Dansk (Danish).....	6
Nederlands (Dutch).....	7
Eesti (Estonian).....	8
Suomi (Finnish).....	9
Français (French) Canada.....	10
Français (French) Européen.....	11
Deutsch (German).....	12
Ελληνικά (Greek).....	13
Magyar (Hungarian).....	14
Italiano (Italian).....	15
日本語 (Japanese).....	16
Polski (Polish).....	17
Português (Portuguese) Brasileiro.....	18
Português (Portuguese) Europeu.....	19
Română (Romanian).....	20
Slovenčina (Slovak).....	21
Español (Spanish) Europeo.....	22
Svenska (Swedish).....	23
Türkçe (Turkish).....	24

English

Intended Use: Embedding media for human or animal tissue.

Storage: Keep Tissue Embedding Media dry, out of direct sunlight and under 86 °F (30°C). Note: Embedding Media are still suitable for use even under adverse conditions where media pellets have fused together. Should this occur, gently tap closed bag on a flat surface to obtain appropriate sized pieces of paraffin..

Hazardous ingredients: Paraplast Plus contains < 2 % DMSO.

Directions For Use: Use in accordance with established tissue embedding and processing procedures.

Melting Point

Paraplast X-tra	52°C
Paraplast, Paraplast Bulk and Paraplast Plus	56°C
Blue Ribbon, Tissue Infiltration Medium, Formula "R", EM-400 Embedding Medium	56 - 57°C
Paraplast High Melt	60°C


Warnings and Precautions:

1. To avoid brittleness and inconsistency in tissue appearance, ensure that wax melt temperature is compatible with tissue.
2. If using separate embedding and infiltration media, ensure compatibility prior to use.

Symbol	Definition
	In Vitro Diagnostic Medical Device
	Catalog number
	Keep Dry
	Keep Away from Sunlight
	Temperature Limitation
	Caution

Symbol	Definition
	Legal Manufacturer
	Authorized Representative in the European Community
	Consult Instructions for Use
	Conformity European
	This Way Up

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Leica Biosystems Richmond, Inc.