Human Uterus

Gill's III Hematoxylin Stain Procedure

0:05

Pint Item#: HXGHE3PT

Liter Item#: HXGHE3LT

Gallon Item#: HXGHE3GAL

5 Gallons Item#: HXGHE35CB

Control Slide(s)	Item#
Trichrome	CST0125P
Skeletal Muscle	CSS0725P
Smooth Muscle	CSS0425P

Required Components (Not Included)	
Differentiating Solution	70% Reagent Alcohol
Bluing Solution	Eosin Y or Aqueous Eosin Y

PRINCIPLE AND RESULTS: This stain is intended for use by laboratory professionals to stain routinely prepared paraffin embedded tissue specimens (in vitro) to identify nuclei and cytoplasm. Nuclei are regressively stained blue to purple and cytoplasm, collagen, and muscle pink to red.

SPECIMEN CRITERIA: Appropriately fixed, paraffin-embedded, 5µm tissue section.

STORAGE AND USAGE NOTES: Store/Use each component according to the temperature and expiration on the label. If using a Xylene Substitute, deparaffinize for 10 minutes each change with agitation and clear for 2 minutes each change with agitation.

PRECAUTIONS: For use by laboratory professionals. See SDS for complete warnings, precautions, hazard and precautionary statements, and disposal information.

STAINING PROCEDURE:

    Color coordinated steps denote stain baths that can be reused during autostainer configuration.

#	Action	With	Heat °C	Time		Details
				Mins	Secs	
1	Deparaffinize	Xylene or Substitute, 2 changes	--	5	--	5 minutes each change or as required if using a xylene substitute.
2	Rinse	Absolute Alcohol, 3 changes	--	1	--	1 minute each change or as required if using graded alcohols.
3	Rinse	95% Reagent Alcohol	--	1	--	1 minute each change with agitation.
4	Rinse	Running Tap Water	--	1	--	
5	Immerse	Gill's III Hematoxylin	--	5	--	Once complete, rinse is running tap water (1 minute) and continue.
6	Immerse	Differentiating Solution	--	--	1	2-3 dips for 1 second each. Once complete, rinse is running tap water (1 minute) and continue.
7	Immerse	Bluing Solution	--	--	30-40	Once complete, rinse is running tap water (2-3 minutes) and continue.
8	Immerse	70% Reagent Alcohol	--	1	--	Agitate.
9	Immerse	Eosin Y or Aqueous Eosin Y	--	--	30-60	Agitate.
10	Dehydrate	Absolute Alcohol, 4 changes	--	1	--	1 minute each change with agitation.
11	Clear	Xylene or Substitute, 3 changes	--	1	--	1 minute each change or as required if using a xylene substitute.
12	Coverslip	Permanent Mounting Media	--	--	--	

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 143-144.
2. With modifications by AMTS R&D Department, 1979-2019.

CE MARKINGS AND DESIGNATIONS:

REF	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
LOT	Batch Code		Use By		Consult Instructions Prior to Use	
 GHS08	Health Hazard					

CONTACT INFORMATION:

American MasterTech Scientific
1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab
2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

1. Sheehan DC Hrapchak BB: Theory and Practice of Histotechnology; 1980, 143-144.
2. With modifications by AMTS R&D Department, 1979-2019.

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