

# ColorBond™

## Tissue Marking Dye Mordant Spray \* 8 OZ. BOTTLE (6 BOTTLES PER CASE)

A specially formulated mordant solution that facilitates marking dye adherence to tissue and helps to prevent bleeding or fading of dyes. Particularly helpful with fatty tissue.

### How does it work?

- ✱ Spray onto dyed tissue as needed for better adherence through processing/staining solutions
- ✱ Facilitates the bonding of marking dyes to tissue
- ✱ Developed to better define margins bleeding, changing color or fading
- ✱ Particularly useful with fatty tissue
- ✱ Especially useful for certain problematic low viscosity dye formulations
- ✱ Packaged in easy to use 8 oz. spray bottles (6 bottles per case)



Cat. **CBM-6** Qty. **6/cs**

### Related Frozen Section Accessories:

- ✱ StatLab Tissue Marking Dyes, 2 oz bottles (7 colors to choose from)  
**Cat. SL662XX-2**
- ✱ Tissue Marking Dye, Acrylic Rack (holds dyes, scalpel blades, applicator sticks, and forceps)  
**Cat. SL667**
- ✱ Tissue Marking Dye Set (includes acrylic tray, marking dyes (7), and applicator sticks)  
**Cat. SL669**
- ✱ Wooden Tip, and Cotton Tipped Applicator Sticks  
**Cat. 826-W/CTA-806W**
- ✱ FormGuard™ Plus, Formalin Neutralizing Pads (8" x 11", 11" x 15", 15" x 20")  
**Cat. FGP0811, FGP1115, FGP1520**
- ✱ Cincinnati Surgical grossing blades, #15, #21, #22, #60  
**Cat. CS-B15T, CS-B21, CS-B22, CS-B60**