



American MasterTech  
scientific laboratory supplies

# 10% Buffered Formalin

## Ready to Use Histological Tissue Fixative

(Rev. 3/8/18)

### 10% Buffered Formalin, The product and its purpose:

10% Buffered Formalin is a Histological tissue fixative intended for In Vitro Diagnostic use. Tissue fixed with 10% Buffered Formalin can be used for routine Special Stains and Immunohistochemical (IHC) stains. Antigen retrieval will be needed for IHC staining. Active ingredients include; Formaldehyde (<5%) and Methanol (<3%).

### 10% Buffered Formalin, Instructions for use:

Use 10% Buffered Formalin as you would any routine histology fixative. Most surgical tissue specimens require 6 to 24 hours for complete fixation. Large tissues may take longer to ensure complete penetration and fixation.

### 10% Buffered Formalin, Storage and Stability:

10% Buffered Formalin should be stored at 2-55 C and is stable under normal laboratory conditions. See product label for shelf life.







### 10% Buffered Formalin, Precautions:

See Safety Data Sheet (SDS) for warnings and precautions.

### 10% Buffered Formalin, Disposal:

Dispose of or neutralize 10% Buffered Formalin in accordance with all local, state, federal, environmental and/or EU regulations.

### CE MARKINGS AND DESIGNATIONS:

<b>REF</b>	Catalogue Number		Temperature Limitation		Manufacturer	American MasterTech Scientific 1330 Thurman St. Lodi, CA 95240 USA Tel 800 860 4073 Fax 209 368 4136
<b>LOT</b>	Batch Code		Use By	<b>EC REP</b>	Representative	Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands
<b>IVD</b>	In Vitro Diagnostic Medical Device		Consult Instructions Prior to Use	<b>CE</b>		
 GHS07	Irritant	 GHS08	Health Hazard			

Copyright © 2015-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.