

# Instructions for Use

## Sterile Screw Top Urine Specimen Containers



CATALOG NUMBER	DESCRIPTION	UNIT OF MEASUREMENT
CT1019-B	4 oz, sterile container	100 per case

### INTENDED USE

Sterile specimen containers are used for the collection and transport of biological specimens, including urine, to ensure there are no outside microorganisms or other contaminants introduced.

### APPLICATIONS

Collection and transport of urine specimens.

### STORAGE AND STABILITY

Storage: Room temperature

### FEATURES

- 1 Made of strong flexible polypropylene material
- 2 Tab type, tamper-evident
- 3 Patient ID label applied to container (also says Sterile and holds the lid closed).

Please contact [tech@statlab.com](mailto:tech@statlab.com) with any additional questions.



#### StatLab Medical Products

2090 Commerce Drive, McKinney, TX 75069

StatLab.com Phone: +1 800-442-3573

