

Instructions for Use

Screw Top StatClick™ Empty Specimen Containers



CATALOG NUMBER	DESCRIPTION	UNIT OF MEASUREMENT
CTL20	20mL, 1 5/8" H x 1 1/4" W	800 per case
CTL40	40mL, 1 9/16" H x 1 11/16" W	400 per case
CTL60	60mL, 2 5/16" H x 1 11/16" W	500 per case
CTL90	90mL, 3" H x 2.3" W	300 per case
CTL120	120mL, 3 1/16" H x 2" W	200 per case

INTENDED USE

StatLab StatClick™ Empty Specimen Containers are designed for safe, leak-resistant transport and storage of fluids. Screw-top containers have a knurled cap and base for easy handling. They are graduated and made of durable clarified polypropylene.

APPLICATIONS

Transport and storage of fluids. StatClick™ containers produce an audible "click" when the leak-proof seal has been properly engaged, assuring accurate closure and maintaining specimen integrity during transport.

STORAGE AND STABILITY

Storage: Room temperature

FEATURES

- 1 Made of durable clarified polypropylene
- 2 StatClick™ container

Please contact tech@statlab.com with any additional questions.



StatLab Medical Products

2090 Commerce Drive, McKinney, TX 75069

StatLab.com Phone: +1 800-442-3573

