

Instructions for Use

BlockVault™ Filing System



CATALOG NUMBER	DESCRIPTION	UNIT OF MEASUREMENT
BVFS-6D	BlockVault™ Cassette Filing System, 6 Drawer	1 each

INTENDED USE

BlockVault™ Filing System is used for the secure, long-term storage of paraffin blocks.

GENERAL INFORMATION

BlockVault™ Filing System offers lightweight, secure, convenient long-term storage of paraffin blocks. The high impact poly construction, six-drawer cabinet accommodates up to 1,500 paraffin blocks. Dimensions are 17 ¼ x 9 ½ x 15 ½.

Please contact tech@statlab.com with any additional questions.



StatLab Medical Products
2090 Commerce Drive, McKinney, TX 75069
StatLab.com Phone: +1 800-442-3573

