

Instructions for Use

Cervical Scraper



CATALOG NUMBER	DESCRIPTION	UNIT OF MEASUREMENT
CT1760	Cervical Scraper, Plastic, 8"	500 per case

INTENDED USE

8" dual end plastic cervical scraper with tapered tip on one end and Ayres tip on the other end for the collection of cell samples.

GENERAL INFORMATION

The dual end plastic cervical scraper can be used to collect specimens from the cervix, or from the tongue, throat, or vulva.

PROCEDURE

Use in accordance with policies and procedures developed by your laboratory.

Please contact tech@statlab.com with any additional questions.



StatLab Medical Products
2090 Commerce Drive, McKinney, TX 75069
StatLab.com Phone: +1 800-442-3573

