

Instructions for Use

Double Edge Grossing Blades



CATALOG NUMBER	DESCRIPTION	UNIT OF MEASUREMENT
DEB-0200	Blades Grossing Double Edged	100 per box
DEB-H	Blade Holder (pencil style) for Double Edged Blades	1 each

INTENDED USE

Double Edged Grossing Blades are made exclusively for histological applications. They are designed for use prior to specimen processing and are used to hand trim specimens. Highly polished and made of carbon steel, they come individually packaged for safer dispensing.

GENERAL INFORMATION

Double Edged Grossing Blades are packaged individually and measure 2 ¼ x ¾. The blades can be securely locked into position with a simple twist of the blade handle.

Double Edge Grossing Blades are best used with catalog number DEB-H, Blade Holder (pencil style) for Double Edged Blades.

DISPOSAL OF BLADE

Blades are for single use and should be disposed of in a sharps disposal container after use.

Please contact tech@statlab.com with any additional questions.



StatLab Medical Products
2090 Commerce Drive, McKinney, TX 75069
StatLab.com Phone: +1 800-442-3573

