

Instructions for Use

Slide Vault Storage Cabinet



CATALOG NUMBER	DESCRIPTION	UNIT OF MEASUREMENT
SVFS-6D	Slide Vault Storage Cabinet	1 each
SVFS-B	Slide Vault Base Unit	1 each

INTENDED USE

Slide Vault Storage Cabinet contains 6 drawers and is used for easy storage and retrieval of slides. Base Unit provides solid support and sits on scratch-proof rubber feet. Each drawer can hold up to 750 slides for a total of up to 4500 slides per cabinet.

GENERAL INFORMATION

Metal, stackable, interlocking, with baked enamel finish. Contains 6 drawers per cabinet with no plastic inserts needed to separate slides. Polished brass knobs facilitate easy storage and retrieval of slides. Up to 10 cabinet units can be safely stacked on a single base unit.

Please contact tech@statlab.com with any additional questions.



StatLab Medical Products
2090 Commerce Drive, McKinney, TX 75069
StatLab.com Phone: +1 800-442-3573

