

MaxShield™ Protective Apparel

High Quality Protection at an Economical Price!



Introducing MaxShield™, a new line of protective apparel made by adhesively bonding a micro porous polyethylene film to a polypropylene spunbonded substrate. The polyethylene film on the outside is extremely durable and provides strength and protection from both solid particulate and liquids. MaxShield is an excellent cost effective alternative to Tyvek products.



MaxShield Apron 2010 Max:

Size:	Qty:	Item#:
--	100	AAA01C

MaxShield Lab Coats w/Elastic Cuffs:

Size:	Qty:	Item#:
Small	30	AAC0102C
Medium	30	AAC0104C
Large	30	AAC0106C
XL	30	AAC0107C
XXL	30	AAC0108C
XXXL	30	AAC0109C

MaxShield Shoe Covers:

Size:	Qty:	Item#:
L-XL	100	AAZ0211C

MaxShield Sleeves:

Size:	Qty:	Item#:
--	100	AAS01C

MaxShield Pants w/Elastic Waist:

Size:	Qty:	Item#:
Small	50	AAP0102C
Medium	50	AAP0104C
Large	50	AAP0106C
XL	50	AAP0107C
XXL	50	AAP0108C
XXXL	50	AAP0109C



American MasterTech
scientific laboratory supplies

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.