

1330 Thurman St. • Lodi, CA. 95240 ☎ 800.860.4073 📠 209.368.4136 🌐 www.americanmastertech.com

Tomahawk[®] Safety Glasses

Providing comfort, protection and style...



American MasterTech
scientific laboratory supplies



Tomahawk[®] safety glasses feature a perfect blend of style and comfort. The lightweight nylon frame and a one-piece polycarbonate lens provide a wide range of vision and angular coverage. Fully adjustable ratchet-action temples provide an excellent fit. Filters out virtually 100% of UV radiation. Meets ANSI Z87.1 standards. Available in blue or black

Color:

Tomahawk[®] Safety Glass Black

Tomahawk[®] Safety Glass Blue

Item#:

AAG03BK

AAG03BL

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.