

Adenna's

Earloop Face Masks



- 4-Fold forms better cone shape for maximum protection
- Fluid Resistant
- Available in Violet, Pink, Blue & White
- Ultra soft inner layer for comfortable wear
- Latex free and fiberglass free
- PVC concealed nose piece minimizes x-ray reflection
- BFE (Bacterial Filtration Efficiency) > 99%
- PFE (Particulate Filtration Efficiency) > 99%

The Adenna 4-Fold Fluid Resistant Earloop Face Masks offer a dependable solution in applications where maximum coverage is preferred. These masks are 100% latex and fiberglass free. The unique and concealed PVC nose piece conforms to your face better and minimizes x-ray reflection. Masks have a very soft inner layer for comfortable wear. PFE (Particulate Filtration Efficiency) > 99% and BFE (Bacterial Filtration Efficiency) > 99%. Ten boxes of 50 masks per case. Available in blue, pink, violet, and white.

Color:	Item #:
Blue/Cs (500)	AAR01BLC
Pink/Cs (500)	AAR01PKC
Violet/Cs (500)	AAR01VIC
White/Cs (500)	AAR01WHC



American MasterTech
scientific laboratory supplies

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.