

encore

by DURAEDGE®

the sharpest, most durable microtome blade ever made.



Discover the performance you demand, with the new Encore low-profile microtome blades by DURAEDGE®. When you compare Encore low-profile blades with blades from other leading competitors, you'll discover how the smallest difference can make a tremendous impact on the quality of your cuts. Encore is the only brand on the market to offer a full 80mm of usable cutting edge. With this full 80mm, you can easily cut up to three standard block widths (24mm) without overlapping. To enhance sharpness, Encore blades are made with an ultra-thin gauge (.010") of steel, which allows the blades to achieve a more acute or steeper angle. Each blade is also uniformly coated with a special proprietary PTFE application process to enhance the blade's durability. The result is longer-lasting, durable blades that can provide consistent, usable sections time after time. The Encore blades are also shown to provide 25-35% more cuts than the leading competitor's blade, which in turn improves productiveness and saves time and money in the long run.



American MasterTech
scientific laboratory supplies

ENCORE LOW-PROFILE MICROTOME BLADES

Pack (50 blades) Case (10 packs)
Item#: BLM0104P Item#: BLM0104C

ENCORE HIGH-PROFILE MICROTOME BLADES

Case (10 packs)
Item#: BLM0106C

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.