



Erado-Sol™

Fast Effective Biological and Chemical Stain Remover

(Rev. 12/6/18)

ERADO-SOL™, PRODUCT INFORMATION:

There is absolutely nothing like it! Erado-Sol rapidly removes most stains from hands, glassware and laboratory equipment. And it's so easy to apply, simply rub in gently and stain is removed in a matter of seconds. Rinse with water. For stubborn (old or set) stains, repeat is necessary.

ERADO-SOL™, REMOVES:

Most medical, histological, and chemical stains including Iodine, Wright -stain, Lugol's stain, and fresh blood.

LABCOATS – HOSPITAL LINENS:

Place a small amount of Erado-Sol on each stain and launder. For stubborn stains, let the liquid soak into stain for 10 to 15 minutes, and then launder.

CONTACT INFORMATION:

American MasterTech Scientific

1330 Thurman St.
Lodi, CA 95240 USA
Tel 800 860 4073
Fax 209 368 4136
www.americanmastertech.com

StatLab

2090 Commerce Drive
McKinney, TX 75069 USA
Tel 800 442 3573
Fax 972 436 1369
www.statlab.com

Copyright ©2011-2019 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.